<u>WNV Pages from the</u> <u>Microbiology Client Services Manual</u>

State of Utah Public Health Laboratory

46 North Medical Drive Salt Lake City, UT 84113-1105

Phone: 801-584-8400 FAX: 801-584-8486

Utah Department of Health

Revised 01/07/05 Page 1

MICROBIOLOGY CLIENT SERVICES MANUAL Utah Public Health Laboratory

GENERAL INSTRUCTIONS

CONTACT US:

ADDRESS, PHONE, FAX, and WEBSITE

Utah Public Health Laboratory 46 North Medical Drive Salt Lake City, UT 84113-1105

Phone: 801-584-8400 FAX: 801-584-8486

Webpage: HTTP://health.utah.gov/els/microbiology

KEY PERSONNEL

Billing

Bob Anderson

Environmental (Water) Microbiology

Sanwat Chaudhuri, Ph.D. -- Section Chief

Microbiology Bureau

Barbara Jepson, MPA, MT(ASCP) -- Bureau Director

Dan Andrews, MS, MT(ASCP) -- Section Chief of Bacteriology,

Food Bacteriology, Mycobacteriology, Parasitology

Norm Brown, BS, MT(ASCP) -- Section Chief of Newborn Screening

Jana Coombs, BS, M/SV (ASCP) - Section Chief of Molecular

Biology, and Bioterrorism Coordinator

Tom Sharpton, MS, SM(ASCP) -- Section Chief of Immunology, Virology

Technical Services

Chris Peper, MT(ASCP) -- Section Chief

REPORTING:

You must supply your correct Customer ID Code to receive test results.

Some mail services and couriers are taking a week or more to get your samples to us. If you are having problems with turn around time for results, check your delivery method. See individual test for specific reporting criteria and methods.

REQUISITIONS:

Blank request forms with your customer ID code are available from Technical Services (also see Appendix B for blank forms WITHOUT the customer ID).

All information must be provided. Incomplete requisitions cannot be processed.

SPECIMEN LABELING: See individual requirements under specific test.

***NOTE: Specimen containers from the State of Utah Public Health Lab have an <u>outdate</u> printed on the label. Do not collect any sample in an outdated container. Call Technical Services at 801-584-8204 for a new container.

We do not supply blood collection tubes.

Revised 10/26/05 Page 2

MICROBIOLOGY CLIENT SERVICES MANUAL Utah Public Health Laboratory

LAB TEST – Immunology or Molecular Biology Sections

TEST West Nile Virus (Human) IgM ELISA

St. Louis Encephalitis Virus (Human) IgM ELISA

METHOD Enzyme Linked Immunosorbent Assay (ELISA)

AVAILABLE Prior to submitting specimen, contact UDOH Epidemiology at

(801)538-6191.

PATIENT PREP Symptoms, vaccinations, and travel history

SPECIMEN Serum or cerebrospinal fluid

COLLECT IN N/A

PROCESSING Serum: refrigerate (freeze if transport delayed)

CSF: refrigerate if transport delayed

TRANSPORT Serum: refrigerate during transport (freeze if transport delayed)

CSF: refrigerate if transport delayed

TIME CRITICAL Within 12 hrs of collection

LABEL Patient's full name or unique ID number, date of collection, and date of

onset of symptoms

REQUISITION Molecular Biology Test Request Form (see form in Appendix B)

TEST COMPLETE 72 hrs after receipt in our lab

RESULTS WNV or SLE antibody detected by ELISA; WNV or SLE not detected by

ELISA

REPORTED Phone, fax, or email, as established with provider

NOTE If initial serum specimen was collected within 9 days of onset of

symptoms, a convalescent serum will be requested for IgM negative tests.

CONTACT (801)584-8449: Jana Coombs or Kim Christensen.

Revised 10/26/05 Page 3

MICROBIOLOGY CLIENT SERVICES MANUAL Utah Public Health Laboratory

LAB TEST – Molecular Biology Section

TEST West Nile Virus, St. Louis Encephalitis Virus, or Western Equine

Encephalitis Virus

METHOD Polymerase Chain Reaction (PCR)

AVAILABLE Contact UDOH Epidemiology at (801)538-6191 or Division of Wildlife

Resources at (801) 538-4767 for submitting avian oral swabs and dead

bird reports.

PATIENT PREP N/A

SPECIMEN Mosquitoes = 10-50 insects, available Mosquito Abatement Districts.

Avian oral swabs.

Bird or horse tissues = 1 cubic centimeter brain, spleen, or heart.

COLLECT IN Mosquitoes = tubes from Mosquito Abatement District.

Swabs = Ziploc bags; outer bag must be clean.

Tissue = sterile, leak proof container.

PROCESSING Keep mosquitoes and tissue samples at 2 - 8 degrees C.

Avian oral swabs at ambient temperature.

TRANSPORT On wet ice or in mailer

TIME CRITICAL Within 48 hrs of collection

LABEL Location and date of collection. Species of source animal.

REQUISITION Molecular Biology Test Request Form (see form in Appendix B)

TEST COMPLETE 48 hrs after receipt in our lab

RESULTS Virus detected by PCR; virus not detected by PCR

REPORTED Mail, e-mail, or fax, as established with provider

NOTE N/A

CONTACT (801)584-8449: Jana Coombs or Kim Christensen.

Revised 10/26/05 Page 4

MOLECULAR BIOLOGY TEST REQUEST FORM

STATE OF UTAH PUBLIC HEALTH LABORATORY

46 NORTH MEDICAL DRIVE SALT LAKE CITY, UTAH 84113-1105

FOR LABORATORY USE ONLY
LAB#:
DATE STAMP:

TELEPHOI	NE: (801) 584-8	400 FAX: (801) 584-	8486				
TESTING WILL <u>NOT</u> BE PERFORMED UNI	ESS SLIP IS <u>CC</u>	MPLETELY FILLED OU	T. PLEASE	PRINT <u>CLE</u> A	RLY FOR ACC	URACY.	
PATIENT INFORMATION:							
Patient Name (Last, First):							
Patient ID #:	DATE OF	DATE OF BIRTH (mm/dd/yy)		SEX:			
	/				М	F	
PROVIDER INFORMATION: Physician:		:		SPECIMEN	N COLLECTIO	N DATE	
Provider Code: Provider F		Provider Phone:			(MM/DD/YY)		
		Email:					
	Secure Fa	x #:		/	/		
[] Serum [] CSF [] Nasopharyngeal swab [] Aspirate (specify):	[] Bordetella pertussis PCR [] West Nile Virus PCR [] St. Louis Encephalitis Virus PCR [] Western Equine Encephalitis PCR [] Human West Nile Virus IgM ELISA [] Human St. Louis Encephalitis Virus IgM ELISA [] SARS PCR [] ORSA PFGE [] Varicella zoster virus PCR [] Vaccinia virus PCR [] Other (specify):						
RESULTS	STATE OF ORIGIN OF PATIENT/SAMPLE						
		ADDITIONAL INFORMATION (List pertinent information including presumptive ID)					

Revised on 1/7/2005 Page 75