

Utah Newborn Screening


SEE BACK FOR BLOOD SPOT COVER

SECOND SCREEN: General Instructions

Collect specimen after 7 days of life.

COLLECTION INSTRUCTIONS

1. Legibly print ALL information in spaces provided using block capital letters.
2. Collect specimen with heel stick. See newborn screening handbook for detailed instructions.
3. Fill all 7 circles.
4. Dry 3-4 hours before mailing.



COLLECT SAMPLE FROM
SHADED AREA.

Note Expiration Date:
Form CANNOT be
used after this date.

Mailing Instructions:

1. When blood is dry, fold stock card (from back of form) over blood spots. The flap should enclose the blood spots and reveal a biohazard symbol.
2. If using the postal service, place form with blood spots covered into envelope.

Mail to: Newborn Screening Laboratory
Utah Department of Health
46 N Medical Dr
Salt Lake City UT 84113-9903

Mother's Information*: This is used to
identify the baby and mother, and to
contact mother if there is a problem.

For more information, call, refer to your
handbook or visit our website:
http://www.health.utah.gov/newbornscreening/HCP_Instructions.htm

Retain this sheet for your records.

XXXXXXX I.D. Number

PEEL AWAY THIS PART HERE

2



FOR UDOH LAB ONLY - DO NOT MARK

XXXXXXX



FOR UDOH LAB ONLY - DO NOT MARK

XXXXXXX

XXXXXXX

XXXXXXX

UTAH DEPARTMENT OF HEALTH
SECOND NEWBORN SCREENING FORM

BLOCK PRINT ALL CAPITALS - COMPLETE ENTIRE FORM

FORM EXPIRES DECEMBER 2007

01 02 2007

Sample collection date MM/DD/YYYY

Medical Record Number

SMITH

Baby's last name

YOUR HOSP

Birthing place/hospital

01 01 2007

Baby's birthdate MM/DD/YYYY

3520

Birthweight (gms)

SMITH

Mother's legal last name

DOE

Mother's maiden name

44 N MEDICAL DR

Mother's mailing address

SLC

City

01 02 1986

Mother's birthdate MM/DD/YYYY

801 584 8256

Mother's Area Code & phone

JIM JONES

Baby's Medical Home: Doctor's Name / Clinic's Name

1234 S MEDICAL DR

Baby's Medical Home: Doctor's / Clinic's Address

SLC

City

UT 84104

State

801 584 8256

Baby's Medical Home: Doctor's / Clinic's Area Code & Phone

RECALL SCREEN MARK ONLY IF INSTRUCTED

Unacceptable 1st 2nd Positive

BELOW FOR UDOH LAB ONLY - DO NOT MARK

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Sample Unacceptable

Fill all 7 circles with blood.
The top two blood spots are
sent off-site and MUST be
filled in addition to the five
below.

FOR UDOH LAB
ONLY: Do not mark or
place labels in any of
these areas.

Sample Collection Date:
Results are based on age of
blood and cannot be
processed without date.

Birthdate: Results are
based on baby's age and
cannot be processed
without birth date.

Birth Weight: Results are
based on birth weight and
cannot be processed
without it.

Medical Home Information:
Who is called in case of an
abnormal screen and where
to send results.

Recall Screen Box:
Mark only if instructed by
program staff. This is used
if another specimen is
needed because of specimen
being unacceptable or
abnormal.

* If baby is adopted, write *adoptive family information*, or adoption agency and contact person. We must be able to identify and connect the information from the second screen with the first screen.