


# Newborn Screening Cards 101:

## Exposing the two additional circles made simple




① Read all instructions on the screening card cover before proceeding. Gently remove cover from screening card

② Remember to check the expiration date of all cards before each use.

③ Fill out card completely, especially the shaded areas.



④ To expose the two additional circles pull on top of protective flap.

**Caution!** Do not separate the two additional circles from the card.



⑤ Make sure to fold the protective flap over the blood spots. Postal regulations require this protective flap to cover the blood spots.

⑥ ATTENTION: The specimens should be sent to the State Lab within 24 hours of collection.



Newborn Screening Follow Up Program  
44 N Medical Drive  
PO Box 144710  
Salt Lake City, UT 84114-4710

For more information  
Call: 801-584-8256

Visit us on the Web:  
[www.health.utah.gov/newbornscreening](http://www.health.utah.gov/newbornscreening)

Email us:  
[newbornscreening@utah.gov](mailto:newbornscreening@utah.gov)