# An Overview of the New Laboratory Information Management System (LIMS) at the Public Health Laboratory

Dan Andrews, MS, MT(ASCP)
Section Chief, Bacteriology/TB Laboratory
Utah State Public Health Laboratory

- Will cover the Bureau of Microbiology only
- History of LIMS in the Bureau
- What previous systems could do
- The New System
- What it does
- What it can do


### **Bureau of Microbiology**

- Bacteriology/Tuberculosis
- Immunology
- Molecular Biology / BT
- Newborn Screening
- Virology

### **Laboratories Using New LIMS**

- Bacteriology/Tuberculosis
- Molecular Biology / BT
- Newborn Screening
- Virology / Immunology
- Toxicology
- Environmental Microbiology

# In the Beginning



#### **Previous LIMS**

- Store patient demographics
- Store patient test results
- Send report to provider
- Send special reports to Public Health Departments\*
- Query database tables

# <87 87 88 89 90 91 92 93 94 95 96 97 98 99 00 01 02 03 04 05 06 07 08</p> Newborn Screening



#### Microbiology



#### LabeWare LIMS

- Store patient demographics
- Store testing data
- Store patient test results
- Store provider information
- Store reagent information
- Store Quality Control information

- Send report to provider
- Send special reports to Public Health Departments
- Quality Control reports
- Query database tables\*

#### **Data screen in AREV**

```
Microbiology Log-
                                                                05/09/2005
 Lab Number:
                                                Date Received:
                  50500450
                                                                         Sex: M
 Name (or ID):
                                                     ATTN: MICROBIOLOGY LAB
                         LDS HOSPITAL
 Provider:
                 LDSH2
 Date Collected: 04/06/2005
                                  DOB: 07/16/1924
                                                Billing Code:
                                                Billing Status:
 Source/Site:
                 ST SPUTUM
   Specifics:
 Test Requested: IM | IDENT, OF MYCOBACTERIA | Acsn by: CR
                                                Data by: CR
                             State:
                                                County:
 Comments:
                 AFB(+)
                                                Verified By:
 Date Reported:
 Tech Code:
 Results:
                 NTBC Acid-fast organism not M. tuberculosis, ID to follo
                        By DNA probe_
 Remarks:
Window Enter
                                  SoftKeys | Relations | MICRO. DAT
                                                                     50500450
|Ed Ins|
                                    Micro.3
                                                                         Level
```

### **Work Card from Bacteriology**

550N

02/02/06 ec'd: 70600217 ab #:

ıme:

Sex: ie: TIMPANOGOS REG. HOSPITAL covider:

ollected: 01/31/06

FECES urce:

ID, BACTERIAL est: R/O SHIGELLA mments:


2/6/04

Mac-> 24HNLF MTO-> -1-H Rapid One

# **Testing Records**



### Labwear LIMS



### **Labware LIMS MIO Detail**

MIO / 1	Version 10	Authorized
Date Completed	08/07/2008 10:36:01 AM	
Date Reviewed	08/08/2008 05:23:47 PM	CKFISHER
Motility result	Negative	
Indol result	Negative	
Ornathine result	Negative	
MIO lot #	HARDY-MIO-08064	
Kovac lot #	HARDY-07354	
Date Received	08/06/2008 02:55:30 PM	
Changed On	08/08/2008 05:24:05 PM	
Date Completed	08/07/2008 10:36:01 AM	

# **MIO QC Detail**

QC_PASS_FAIL / 1	
Lot QC result	Pass
QC organism 1 ID	Escherichia coli ATCC 25922
QC organism 2 ID	Klebsiella pneumoniae ATCC 33495
QC organism 1 result	motile+, indole-red, ornithine-purple
QC organism 2 result	non-motile, indole-yellow, ornithine-yellow
Sterility Test result	No growth seen



Utah Public Health Laboratories 45 North Medical Drive Salt Lake City, Utah 84113-1105 Telephone: (801) 584-8400 Fax: (801) 584-8501

PROVIDER CODE: IMC2

INTERMOUNTAIN CENTRAL LAB-MICRO

MICROBIOLOGY LAB

PO BOX 57970

MURRAY, UT 84157-0970

Patient Informatio

Patient Name:

Sex: FEMALE

Date Of Birth:

Age: 78 State of Origin: UT Order#: 0000117170

Sample Collection Information

Lab Number: 148423

Source: Stool Specifics:

Date Collected: 07/18/2008 Date Received: 08/06/2008 Date Released: 08/08/2008

Date Reported: 08/09/2008

Test Results

Test Run: Baoterial ID/Reference Collection

Result: SHIGELLA connel

INTERMOUNTAIN CENTRAL LAB-MICRO MICROBIOLOGY LAB PO BOX 57970 MURRAY, UT 84157-0970

### **Pending Report**

8/7/2008

#### DAILY BACTERIOLOGY PENDING REPORT

Login Dates From: 06-Aug-2008 To: 06-Aug-2008

Samp Numb	Patient Name	<u>Sex</u>	Age	Cust ID and Desc	Collected	Source	Comments/Presumptive
148158		F	5	LRMC: LOGAN REG. HOSPITAL	08/02/2008	NOT_SPECIF	CAMPYLOBAC
148420		F	3	ARUP1: ARUP MICROBIOLOGY	08/04/2008	NOT_SPECIF	SHIGELLA
148421		M	53	LRMC: LOGAN REG. HOSPITAL	07/26/2008	NOT_SPECIF	CAMPYLOBAC
148422	<del>-</del>	M	0	IMC2: INTERMOUNTAIN CENTRAL LAB-MICRO	08/01/2008	STOOL	SALMONELLA
148423	<del>-</del>	F	78	IMC2: INTERMOUNTAIN CENTRAL LAB-MICRO	07/18/2008	STOOL STOOL	SHIGELLA
148481		F	66	SMH1: ST MARK'S HOSPITAL	08/01/2008	STOOL	SALMONELLA

# **Daily Positive Report**

8/9/2008 DAILY BACTERIOLOGY POSITIVE REPORT								
		Login Dates From:		ogin Dates From:	08-Aug-2008	To: 08-Aug-2008		
Samp Numb	Patient Name	Sex	Age	<u>DOB</u>	Cust ID and Desc	Collected	Source	Result Description
142497		M	24	07/13/1984	QUES: QUEST DIAGNOSTICS INC.	07/22/2008	STOOL	SALMONELLA enteritidis
144019		F	35		IMC2: INTERMOUNTAIN CENTRAL LAB-MICRO	07/27/2008	URINE	SALMONELLA saintpaul
147066	-	F	83	02/09/1925	ARUP1: ARUP MICROBIOLOGY	07/23/2008	STOOL	SHIGELLA flexneri
147900		F	14	09/22/1993	PCMC1: PRIMARY CHILDREN'S MEDICAL CENTER	08/03/2008	CERVIX	NEISSERIA gonorrhoeae
148423		F	78	12/20/1929	IMC2: INTERMOUNTAIN CENTRAL LAB-MICRO	07/18/2008	STOOL	SHIGELLA sonnei
149525	-	М	40	11/25/1967	ARUP1: ARUP MICROBIOLOGY	07/27/2008	SPUTUM	LEGIONELLA pneumophila group 1
149738		М	2	02/13/2006	ACLL: LABORATORY CORP. OF AMERICA	08/01/2008	STOOL	ESCHERICHIA coli, O157:H7

### Ad hoc Queries

- In Labware LIMS all queries must be written by a programmer.
- We have a selection of basic queries that we can work with.
- If a request for information does not fit into one of our choices we have to request programming time to have the query written.

#### What's Next

- Direct transfer of data into provider and health department database (HL7).
- Providers will be able to order tests electronically.