



## ***Chemical and Environmental Laboratory***



# **Client Services Manual**

Revision 1 March 2005  
Revision 2 June 2013

**Utah Public Health Laboratory**

## State of Utah Public Health Laboratory

## CONTACT US:

**Phone:** (801) 965-2400


**FAX:** (801) 584-8251

**Website:** <http://health.utah.gov/lab>

After hours, Holidays, and Emergency: 1-888-374-8824

For Scheduling Sample Testing: 801-965-2405 or 801-965-2470

**MAP OF LOCATION:**



Bret Van Ausdall, PhD: Laboratory Chief

## David Dick: Sample Receiving

Steve Dickson: Manager of Inorganic Chemistry section

Kyle Ashby, PhD: Manager of Organic Chemistry, Radiochemistry, and Water Microbiology Section

Robert Lo: Manager of Metals section.

Manual Edited by Jenna Moffatt, June 2013

## **USE OF THIS CLIENT MANUAL**

This manual is designed to provide a resource for sample collection and available testing information. Its contents are not to be used for regulatory purposes other than providing proper sample collection and preservation information. If you have any questions, please contact Bret Van Ausdall at (801) 965-2458 or (801) 965-2400.

For your benefit, there is a general description and list of tests for the Clean Water Act (CWA), Safe Drinking Water Act (SDWA), and the Resource Conservation Recovery Act (RCRA) on page 7.

To find a specific test, there is a Table of Contents on pages 4-5 and an Index on pages 76-77. There is also a Program Specific Index for each of the above mentioned acts that list each specific test covered by the act. These are located on pages 47-70. A special tests section for Unregulated Contaminant Monitoring Rule 3 (UCMR3) is located on pages 46-47. The Table of Contents is organized alphabetically by contaminant within each analytical section. Follow the page number to locate more details about each contaminant. The index in the back of the manual lists each test alphabetically along with their method numbers. Program specific indices are organized by the specific program, giving the analytes and tests codes.

Once a contaminant has been selected, specific information and testing instructions are located on the page dedicated to the test. The section (organic, inorganic, etc.) is listed at the top of the page. Next, the name and test codes are listed. Analytes being tested for and, if relevant, the grouping of the chemical will follow. To finish the listing, there are the instructions for collection, required types and number of bottles, preservative and handling information, and method technology. Utilize this information to gain a greater understanding of the specific test that is associated with each contaminant.

**All information should be checked and confirmed to be accurate before testing. Instructions should be followed as stated and samples returned to the lab for processing.**

# ***Chemical and Environmental Laboratory – Client Services Manual***

State of Utah Public Health Laboratory

## **TABLE OF CONTENTS**

<u>Analytical Services</u>	<u>Page</u>
Organic Chemistry .....	
BTEX N.....	13
Carbamates .....	11
Glycols.....	15
Haloacetic Acids (HAA).....	8
Herbicides.....	11
Oil and Grease.....	14
Pesticides and SVOCs.....	10
Pesticides (Phase II / Phase V).....	12
Polychlorinated Biphenyls and Organochlorine Pesticides .....	14
Surfactants.....	15
Total Petroleum Hydrocarbons (TPH) .....	13
Toxic Characteristic Leaching Procedure (TCLP-Organics).....	16
Trihalomethanes (THM) .....	8
Volatile Organic Compounds (VOCs).....	9
Inorganic and Metal Chemistry.....	
Alkalinity .....	30
Ammonia .....	22
Annual Inorganics and Metals (18 parameters) .....	18
Arsenic.....	19
Biological Oxygen Demand (BOD) .....	32
Bromate .....	28
Bromide .....	28
Chemical Oxygen Demand (COD) .....	33
Chemistry Types (analytical groupings).....	35
Chlorate .....	28
Chloride .....	28
Chlorite .....	28
Chlorophyll-A .....	33
Chromium-VI.....	36
Color.....	27
Conductivity .....	26
Corrosivity.....	26
Cyanide .....	29
Fluoride.....	30
Hardness .....	27
Lead and Copper .....	18
Mercury.....	19
Metals, Total and Dissolved (analytical groupings).....	20
Digestion Methods for Metals (non-drinking water).....	21
New Drinking Water Source (46 parameters) .....	17
Nitrate and Nitrite .....	23
Nitrite (only) .....	23
Nitrogen, Total Kjeldahl (TKN) .....	25
Nitrogen, Total Dissolved (DTN).....	25
Nutrients, Total and Dissolved (analytical groupings) .....	24
Odor (Threshold Odor Number, TON).....	27
Perchlorate .....	29
Periphyton.....	33

# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## TABLE OF CONTENTS

Inorganic and Metal Chemistry (continued)	
pH.....	25
Phosphate .....	22
Silica.....	34
Sulfate .....	23
Sulfide.....	30
Settable Solids (SS).....	34
Total Dissolved Solids (TDS) .....	34
Total Organic Carbon (TOC).....	31
Total Suspended Solids (TSS).....	34
Total Volatile Solids (TVS) .....	34
Toxic Characteristic Leaching Procedure (TCLP-Metals) .....	37
Turbidity.....	32
UV254.....	31
Environmental Microbiology .....	
<i>Cryptosporidium</i> and <i>Giardia</i> .....	39
Heterotrophic Plate Count.....	38
<i>Enterococcus</i> .....	40
<i>Legionella</i> .....	40
Total Coliform, Fecal Coliform and <i>E. coli</i> .....	38
Radiologic Chemistry .....	
Gross Alpha.....	41
Gross Beta.....	41
Gross Gamma .....	43
Radium-226 .....	42
Radium-228 .....	42
Radon.....	43
Uranium.....	36
Specialty Analyses .....	
Geosmin and MIB .....	44
Lead (air samples) .....	44
UCMR3.....	45
Program Specific Indexes .....	
Clean Water Act (CWA) .....	47
Clean Water Act (CWA) – EPA 608-PCB/OcPest.....	49
Clean Water Act (CWA) – EPA 624-VOC.....	50
Clean Water Act (CWA) – EPA 625-SVOC .....	53
Resource Conservation and Recovery Act (RCRA) – 8081 .....	56
Resource Conservation and Recovery Act (RCRA) – 8151 .....	57
Resource Conservation and Recovery Act (RCRA) – 8260 .....	58
Resource Conservation and Recovery Act (RCRA) – 8270 .....	60
Safe Drinking Water Act (SDWA) .....	63
Safe Drinking Water Act (SDWA) – EPA 524.2-VOC.....	66
Safe Drinking Water Act (SDWA) – EPA 525.2-SVOC .....	69
Appendix-Additional Tests .....	
Selenium Hydride.....	71
Test Request Forms.....	
Instructions.....	72
Index.....	75
General Request Forms .....	77
Chain of Custody.....	80

## ***Chemical and Environmental Laboratory – Client Services Manual***

State of Utah Public Health Laboratory

### **GENERAL LAB PRACTICES AND POLICES:**

Our laboratory is responsible for the receipt and testing of chemicals and environmental microbial contaminants in drinking water, wastewater, environmental soils and hazardous waste. Testing is done to ensure compliance with health and safety standards established by Federal and Utah State agencies. Services are provided to Utah DEQ, public water and wastewater utilities, local health departments, as well as other state and federal agencies.

Laboratory business hours are Monday-Friday 8:00 AM- 5:00 PM. We are closed on holidays with the exception of emergencies. Our staff are on-call 24/7 to handle emergency situations, please call 1-888-374-8824 in case of emergency.

Labeled sample collection materials such as water bottles, glass vials, and solids containers may be obtained from the laboratory (free of charge). The laboratory tests each lot of containers to ensure they are free of contamination. The tested lot number appears on each container label.

A complete test sample submission form (provided when containers are picked up and also included at the end of this manual) should accompany all samples. Tests, which are needed by the client but are not performed in our laboratory, will be subcontracted to a commercial laboratory with approval from the client.

Test results are provided to the submitting client, to other individuals as authorized by the submitting client, and to state and/or federal regulatory agencies as required by law. Fees for laboratory services will be charged to the submitting client.

## ***Chemical and Environmental Laboratory – Client Services Manual***

State of Utah Public Health Laboratory

### **Federal Act Summary**

The Chemical and Environmental Laboratory provides testing to supply the Clean Water Act (CWA), Resource Conservation and Recovery Act (RCRA), and the Safe Drinking Water Act (SDWA).

#### **Clean Water Act**

The Clean Water Act (CWA) is a federal law, which regulates discharges into the navigable waters of the nation. The Utah Division of Water Quality (DWQ) within the Department of Environmental Quality's mission is to protect, maintain and enhance the quality of Utah's surface and underground waters for appropriate beneficial uses; and to protect the public health through eliminating and preventing water related health hazards which can occur as a result of improper disposal of human, animal, or industrial wastes while giving reasonable consideration to the economic impact.

The Utah Public Health Laboratory cooperates with the Division and performs requested water and soil analyses using EPA approved methods.

#### **Resource Conservation and Recovery Act (RCRA) Summary Hazardous Waste and Emergency Response**

The Resource Conservation and Recovery Act (RCRA) is a federal law which defines solid and hazardous wastes and defines requirements for their proper disposal. Both the Utah Division of Solid and Hazardous Waste (DSHW) and the Utah Division of Environmental Response and Remediation within the Department of Environmental Quality protect public health and the environment by ensuring proper management of solid and hazardous wastes and through cleanup of chemically contaminated sites and proper management of underground storage tanks.

The Utah Public Health Laboratories use methods outlined in EPA SW-846 or from other approved sources to test for RCRA contaminants.

#### **Safe Drinking Water Act Summary**

The Safe Drinking Water Act (SDWA) is a federal law which applies to all public water systems. Water systems on tribal lands are regulated by EPA Region 8 out of Denver. All other public water systems are regulated by the Utah Division of Drinking Water (DDW) within the Utah Department of Environmental Quality. DDW's mission is to protect the public against waterborne health risks through education, assistance, and oversight.

The Unified State Laboratories: Public Health is approved to conduct testing of water collected from public water systems according to EPA and State rules. It is not necessary to know the names of each analyte to be requested, but general categories of analytical testing are listed below and can be selected on analysis request forms (see last pages of this manual).

Our lab is endorsed by the EPA to carry out tests for the contaminants that are outlined by these three federal regulations. In the back index, you can look up a contaminant by its specific act and then alphabetically by analyte. There are also method numbers for each analyte and a page.

## ORGANIC CHEMISTRY

Name: **Haloacetic Acids (HAA)**

Test Code: **SM 6251B-HAA**



Grouping: **Disinfection By-Products**

Application: Water systems using chlorine or bromine for disinfection.

Analytes: Dibromoacetic acid  
Dichloroacetic acid  
Monobromoacetic acid  
Monochloroacetic acid  
Trichloroacetic acid

Instructions for Collection:

Allow sample tap to run for a few minutes until water temperature stabilizes. Slowly fill vials to top of container, be sure not to over fill to prevent loss of preservative. Check for air bubbles by inverting, and fill the remaining portion if bubbles found. There must be no headspace.

Required Containers/Volume: 3/40 mL vials

Preservative and Handling: 65 mg  $\text{NH}_4\text{Cl}$ , refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, holding times, 14 days to extract and 7 days to analyze

Method Technology: Liquid-liquid extraction followed by GC detection

Name: **Trihalomethanes (THM)**

Test Code: **EPA 524.2-THM**



Grouping: **Disinfection By-Products**

Application: Water systems using chlorine or bromine for disinfection.

Analytes: Bromodichloromethane  
Bromoform  
Chlorodibromomethane  
Chloroform

Instructions for Collection:

Allow sample tap to run for a few minutes until water temperature stabilizes. Slowly fill vials to top of container, be sure not to over fill to prevent loss of preservative. Check for air bubbles by inverting, and fill the remaining portion if bubbles found. There must be no headspace.

Required Containers/Volume: 3/40 mL vials

Preservative and Handling: 4 mg  $\text{Na}_2\text{S}_2\text{O}_3$ , refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 14 days.

Method Technology: Purge and trap technique followed by GCMS detection





## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### ORGANIC CHEMISTRY

Name: **Volatile Organic Compounds (VOCs)**

Test Codes: **EPA 524.2-VOC, EPA 624-VOC, EPA 8260-VOC or EPA 1666- PMI-VOC**



For lists of analytes see Program Specific Indexes

Clean Water Act (CWA).....page 49

Resource Conservation and Recovery Act (RCRA).....page 56

Safe Drinking Water Act (SWDA).....page 66

#### NOTE:

**Prior to collection of EPA 624, EPA 8260, or EPA 1666 samples, analysis must be scheduled at (801) 965-2400.**

**If using Encore sampler, please bring to lab as soon as possible (within 48 hours).**

#### Instructions for Collection:

**Caution: For surface waters and groundwater samples, check to see if HCl preservative reacts with source water (foams, effervesces, etc.) If reacts, do not add HCl.** For aqueous samples, fill vials with sample source to top of container. Add 2 drops HCl acid to each vial, more drops needed if highly buffered source. Be sure not to over fill to prevent loss of preservative. Check for air bubbles by inverting, and fill the remaining portion if bubbles found. There must be no headspace.

#### Required Containers/Volume:

Drinking Water – 3/40mL vials, plus 1 trip blank  
Surface Water – 4/40 mL, plus 1 trip blank  
\*must use 4/40 mL vials—special prepared for chlorinated sites  
Groundwater – 4/40 mL vials, plus 1 trip blank, prepared and delivered with sample bottles.  
\*must use 4/40 mL vials—special prepared for chlorinated sites  
Soil, solids, sludge – 4 oz. glass container with Teflon-lined lid

#### Preservative and Handling:

**\*24 hour holding for water samples if there is no preservative added**

Drinking Water – 25 mg ascorbic acid, HCl to pH <2  
Surface Water and Method 1666– 10 mg Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> for chlorinated sites and HCl if needed  
Groundwater – 10 mg Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> for chlorinated sites and HCl if needed  
All sample types, refrigerate or store on ice and do not allow to freeze, 4-6°C  
**recommend next day receipt at lab**, total holding time 14 days, except Soil 14 days to extract and 14 days to analyze

#### Method Technology:

Purge and trap technique followed by GCMS detection

#### Preparation Method:


Water (EPA 624 and EPA 8260): EPA 3050  
Soil (EPA 8260): EPA 3585

## Chemical and Environmental Laboratory – Client Services Manual


State of Utah Public Health Laboratory

### ORGANIC CHEMISTRY

Name: **Pesticides and other Semi-Volatile Organic Compounds (SVOCs)**  
Test Codes: **EPA 525.2-SVOC** , **EPA 625-SVOC** or **EPA 8270-SVOC**



EPA 525.2-SVOC  
(Drinking Water)



EPA 625-SVOC  
(Surface Water)



EPA 8270-SVOC  
(Groundwater and  
Solid Samples)

For lists of analytes see Program Specific Indexes

Clean Water Act  
(CWA).....page 53

Resource Conservation and  
Recovery Act  
(RCRA).....page 60

Safe Drinking Water Act  
(SWDA).....page 69

#### NOTE:

**Prior to collection of EPA 525.2, EPA 625 or EPA 8270 samples, analysis must be scheduled at (801) 965-2400.**

#### Instructions for Collection:

For drinking water samples, pour the small vial of acid into each sample bottle. Allow source to flow for a few minutes until water temperature stabilizes, **do not use Tygon tubing**. Slowly fill bottles to top of container, be sure not to over fill to prevent loss of preservative.

#### Required Containers/Volume:

Drinking Water – 2/1L amber glass bottles + acid vials  
Surface Water – 2/1L amber glass bottles  
Groundwater – 2/1L amber glass bottles  
Soil, solids, sludges – 4 oz. amber glass container with Teflon-lined lid

#### Preservative and Handling:

Drinking Water – HCl to pH <2, 50 mg sodium sulfite  
Surface Water – No preservative  
Groundwater and Solids – No preservative  
All sample types, refrigerate or store on ice 4-6°C and do not allow to freeze, **recommend next day receipt at lab**, holding times:  
EPA 525.2 -14 days to extract and 30 days to analyze  
EPA 625 and EPA 8270 – 7 days to extract and 40 days to analyze  
Solids - 14 days to extract and 40 days to analyze

#### Method Technology:

Liquid-solid extraction followed by GCMS detection

#### Preparation Method:

Groundwater (EPA 8270)-EPA 3510  
Soil-(EPA 8270)-EPA 3550



## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### ORGANIC CHEMISTRY

Name: **Herbicides (Chlorinated Organic Acids)**

Test Code: **EPA 515.1-Herbicides or EPA 8151-Herbicides**



Analytes: 2,4-D  
2,4,5-TP (Silvex)  
Dalapon  
Dicamba  
Dinoseb  
Pentachlorophenol  
Picloram

Instructions for Collection:

If applicable, allow sample tap to run for a few minutes until water temperature stabilizes. Slowly fill bottles to top of container, be sure not to over fill to prevent loss of preservative. Collect per project sampling plan.

Required Containers/Volume:

Water – 2/1L amber glass bottles  
Soils, Solids, and Sludge– 4 oz. amber glass container with Teflon-lined lid

Preservative and Handling:

30 mg  $\text{Na}_2\text{S}_2\text{O}_3$  if chlorinated, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, holding times, 14 days to extract and 28 days to analyze

Method Technology:


EPA 515.1 and 8151 (Water)- Liquid-liquid extraction followed by GC detection EPA 8151-(Soils) Solid-liquid extraction, then liquid-liquid extraction followed by GC detection

Prep Method:

EPA 515.1

Name: **Carbamates (Insecticides and Aldicarbs)**

Test Code: **EPA 531.1**



Analytes: 3 - Hydroxycarbofuran  
Aldicarb (Temik)  
Aldicarb sulfone  
Aldicarb sulfoxide  
Carbaryl (Sevin)  
Carbofuran (Furadan)  
Methomyl  
Oxamyl (Vydate)

Instructions for Collection:

Allow sample tap to run for a few minutes until water temperature stabilizes. Slowly fill vial to top of container, be sure not to over fill to prevent loss of preservative. Check for air bubbles by inverting, and fill the remaining portion if bubbles found. There must be no headspace.

Required Containers/Volume:

40 mL amber glass vial

Preservative and Handling:

1.2 mL Monochloroacetic acid, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Method Technology:

Aqueous injection HPLC with post column derivation

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### ORGANIC CHEMISTRY

Name: **Pesticides (Phase II / Phase V)**

Test Codes: **VOCs\* , Pesticides , SVOCs , Herbicides , and Carbamates**



Grouping: **Organic Testing**

Application: Drinking Water

For lists of analytes see Program Specific Indexes

Safe Drinking Water Act (SWDA).....page 63

**NOTE:** Prior to collection of **Pesticides (Phase II / Phase V) Organic Testing**, analysis must be scheduled at (801) 965-2400.

\* **Pesticides (Phase II / Phase V)** may be performed with or without VOC testing

#### Instructions for Collection:

For VOC vials, add 2 drops HCl acid to each vial, more drops needed if highly buffered source. For EPA 525.2-SVOC, pour entire small vial of acid into each sample bottle.

Allow sample tap to run for a few minutes until water temperature stabilizes. Slowly fill to top of containers, be sure not to over fill to prevent loss of preservatives. Check for air bubbles by inverting, and fill the remaining portion if bubbles found. There must be no headspace.

#### Required Containers/Volume:

VOCs - 3/40mL vials, plus 1 trip blank  
Pesticides and SVOCs - 2/1L amber glass bottles + acid vials  
Herbicides - 2/1L amber glass bottles  
Carbamates - 40 mL amber glass vial

#### Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze,  
**recommend next day receipt at lab**, total holding times:  
VOCs - 14 days  
Pesticides and SVOCs - 14 days to extract and 30 days to analyze  
Herbicides - 14 days to extract and 28 days to analyze  
Carbamates - 28 days

#### Methods Technologies:

VOCs - Purge and trap technique followed by GCMS detection  
SVOCs - Liquid-solid extraction followed by GCMS detection  
Herbicides - Liquid-liquid extraction followed by GC detection  
Carbamates - Aqueous injection HPLC with post column derivation

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### ORGANIC CHEMISTRY

Name: **Total Petroleum Hydrocarbons (TPH)**

Test Code: **EPA 8015-TPH**



Analytes: GRO (C6-C10) – Gasoline Range Organics  
DRO (C10-C28) – Diesel Range Organics  
ORO (C28-C35) – Oil Range Organics

-----  
Total TPH (GRO + DRO + ORO)

**Instructions for Collection:**

If applicable, allow sample tap to run for a few minutes until water temperature stabilizes. Slowly fill vials to top of container. Check for air bubbles by inverting, and fill the remaining portion if bubbles found. There must be no headspace.

**Required Containers/Volume:**

Water – 3/40 mL vials

Soil – 4 oz. glass container with Teflon-lined lid

**Preservative and Handling:**

No preservative, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, holding times, 14 days to extract and 40 days to analyze

**Method Technology:**

Extraction technique (liquid-liquid, waste dilution) followed by GC/FID detection.

**Preparation Method:**

Water-EPA 3510; Soil- EPA 3580

Name: **Benzene, Toluene, Ethylbenzene, Xylene, and Naphtalene (BTEX N)**

Test Code: **SM 8260B BTEX N**



Analytes: Benzene  
Toluene  
Ethylbenzene  
Xylene  
Naphtalene

**Note: If using Encore sampler, please bring to lab as soon as possible (within 48 hours).**

**Instructions for Collection:**

If applicable, allow sample tap to run for a few minutes until water temperature stabilizes. Slowly fill vials to top of container. Check for air bubbles by inverting, and fill the remaining portion if bubbles found. There must be no headspace.

**Required Containers/Volume:**

Water – 3/40 mL vials

Soil – 4 oz. glass container with Teflon-lined lid

**Preservative and Handling:**

No preservative, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 14 days

**Method Technology:**

Purge and trap technique followed by GC/MS detection

**Preparation Method:**

Water-EPA 5030; Soil-EPA 3585



## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### ORGANIC CHEMISTRY

Name: **Oil and Grease**

Test Code: **EPA 1664-O/G**



Analyte: Total O/G (Oil and Grease)

**Will be contracted to an outside laboratory with client approval.**

Instructions for Collection:

If applicable, allow source to flow for a few minutes until water temperature stabilizes. Slowly fill containers to top, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

Water – 1 L glass container with Teflon-lined lid

Soil – 4 oz. glass container with Teflon-lined lid

Preservative and Handling:


H<sub>2</sub>SO<sub>4</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days from collection.

Method Technology:

Liquid-solid extraction followed by Gravimetric detection

Name: **Polychlorinated Biphenyls (PCBs) and Organochlorine Pesticides (OcPest)**

Test Codes: **EPA 608-PCB/OcPest or EPA 8081-OcPest and 8082-PCBs**



For lists of analytes see Program Specific Indexes

Clean Water Act (CWA).....page 49

Resource Conservation and Recovery Act (RCRA).....page 56

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

Surface Water – 1 L amber glass bottle

Groundwater – 1 L amber glass bottle

Soils, solids, sludge – 4 oz. amber glass container with Teflon-lined lid

Oil – 4 oz. amber glass container with Teflon-lined lid

Preservative and Handling:

No preservative, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, holding times, 7 days to extract and 40 days to analyze, except Soil 14 days to extract and 40 days to analyze

Method Technology:

Extraction technique (Liquid-liquid, sonication, waste dilution) followed by GCMS detection

Preparation Method:

Water- EPA 3510; Oils-EPA 3580; Soil-EPA 3550

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### ORGANIC CHEMISTRY

Name: **Glycols (Ethylene and Propylene Glycols)**

Test Code: **EPG**



Analytes: Ethylene Glycol  
Propylene Glycol

**Will be contracted to an outside laboratory with client approval. Bottles are to be obtained from contracted lab.**

Instructions for Collection:

If applicable, allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

Water – 2/40 mL vials

Soil – 4 oz. glass container with Teflon-lined lid

Preservative and Handling:


No preservative, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Method Technology:

Aqueous injection with GCFID detection

Name: **Surfactants**

Test Code: **SM 5540C-Surfactants**



Grouping: **New Drinking Water Source**

Analyte: Total Surfactants

**NOTE:** Prior to collection of Surfactants Testing, analysis must be scheduled at (801) 965-2400.

Instructions for Collection:

Allow sample tap to run for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

1 L amber glass bottle

Preservative and Handling:

No preservative, refrigerate or store on ice and do not allow to freeze, **recommend same day receipt at lab**, total holding time 48 hours

Method Technology:

Extraction technique followed by MBAS detection


# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## ORGANIC CHEMISTRY

Name: **RCRA Toxic Characteristic Leaching Procedure (TCLP)**

Test Codes: **TCLP-Organics (VOCs, SVOCs, Pesticides, and Herbicides)**



Grouping: **Organic Testing (RCRA)**

**NOTE:** Prior to collection of EPA 1311-TCLP (water and soil) or EPA 1332-TCLP (oil), analysis must be scheduled at (801) 965-2400.

### Analytes:

#### Volatile Organic Compounds (VOCs)

D029 1,1 - Dichloroethene  
D028 1,2 - Dichloroethane  
D027 1,4 - Dichlorobenzene  
D018 Benzene  
D019 Carbon Tetrachloride  
D021 Chlorobenzene  
D022 Chloroform  
D035 Methyl Ethyl Ketone  
D039 Tetrachloroethene  
D040 Trichloroethene  
D043 Vinyl Chloride

#### Pesticides

Chlordane  
Endrin  
Heptachlor  
Heptachlor Epoxide  
Lindane  
Methoxychlor  
Toxaphene

#### Herbicides

2,4 - D  
2,4,5 - TP (Silvex)

#### Semi-Volatile Organic Compounds (SVOCs)

D030 2,4 - Dinitrotoulene  
D041 2,4,5 - Trichlorophenol  
D042 2,4,6 - Trichlorophenol  
D032 Hexachlorobenzene  
D033 Hexachlorobutadiene  
D034 Hexachloroethane

D024 m - Cresol  
D036 Nitrobenzene  
D023 o - Cresol  
D025 p - Cresol  
D037 Pentachlorophenol  
D038 Pyridine

### Instructions for Collection:

Slowly fill to top of container.

### Required Containers/Volume:

Liquids – 4 L amber glass bottle  
Soil – 4 oz. amber glass container with Teflon-lined lid  
Oil – 4 oz. amber glass container with Teflon-lined lid

### Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, holding times:  
SVOCs - 7 days to TCLP and 40 days to analyze  
VOCs - 14 days to TCLP and 40 days to analyze  
Pesticides - 7 days to extract and 40 days to analyze  
Solids - 14 days to extract and 40 days to analyze  
Herbicides - 14 days to extract and 28 days to analyze

### Method Technology:

Leaching procedure and extraction technique followed by GCMS detection



# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## INORGANIC AND METAL CHEMISTRY

Name: **New Drinking Water Source or “New Source Inorganics and Metals”**  
(46 parameters)

Test Code: **Type 7 (Total Inorganics and Metals Chemistry)**



D = Dissolved

T = Total

L = Lab measured

F = Field measured

Analytes:	Odor (TON)	D-Calcium
	<u>Surfactant</u>	D-Magnesium
	Alkalinity	D-Potassium
	Bicarbonate	<u>D-Sodium</u>
	Carbon Dioxide	T-Aluminum
	Carbonate	T-Antimony
	Carbonate (CO <sub>3</sub> ) Solids	T-Arsenic
	Chloride	T-Barium
	Corrosivity	T-Beryllium
	Fluoride	T-Boron
	Hardness	T-Cadmium
	Hydroxide	T-Chromium
	L-pH or F-pH	T-Copper
	L-Specific Conductivity	T-Iron
	Silica	T-Lead
	Total Dissolved Solids (TDS)	T-Manganese
	Total Suspended Solids (TSS)	T-Mercury
	Turbidity (NTU)	T-Nickel
	<u>Cyanide</u>	T-Selenium
	Sulfate	T-Silver
	Ammonia	T-Thallium
	NO <sub>2</sub> +NO <sub>3</sub>	<u>T-Zinc</u>
	T-Phosphate	Total Organic Carbon (TOC)
	<u>Color</u>	
Total Chemistry	1 L plastic	
Nutrient Analysis	500 mL plastic	
		Filtered Metals
		250 mL plastic
		Total Metals
		250 mL plastic

**NOTE:** Prior to collection of New Drinking Water Source samples, analysis must be scheduled with the laboratory at (801) 965-2400.

Instructions for Collection: Allow sample tap to run for a few minutes until water temperature stabilizes. Slowly fill to top of bottles, be sure not to over fill to prevent loss of preservatives.

Required Containers/Volume:

- 1 L Odor – amber glass bottle
- 1 L Surfactant – amber glass bottle
- 1 L Total Chemistry – unpreserved plastic bottle
- 500 mL Cyanide – plastic bottle
- 500 mL Nutrient Analysis – plastic bottle
- 250 mL Color – plastic bottle
- \*250 mL Filtered Metals – plastic bottle
- \*250 mL Total Metals – plastic bottle
- 200 mL TOC – amber glass bottle

**\*\*Metals Preserve with HNO<sub>3</sub> to pH <2.**

Preservative and Handling: Refrigerate or store on ice and do not allow to freeze, **recommend same day receipt at lab, must be received at lab within 24 hours of collection time**

Methods Technologies: Electrometric, Nephelometric, Gravimetric, Titration, Flow Injection Colorimetry, Ion Chromatography, ICP, and ICPMS

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Annual Inorganics and Metals (18 parameters)**

Test Code: **Type 9 (Primary Inorganics and Metals Chemistry)**



Type 9 (Primary Inorganics and Metals Chemistry)

Analytes:	Cyanide	Fluoride
	Turbidity (NTU)	Total Dissolved Solids (TDS)
	T-Arsenic	T-Barium
	T-Beryllium	T-Cadmium
	T-Chromium	T-Copper
	T-Lead	T-Mercury
	T-Nickel	T-Selenium
	T-Antimony	T-Thallium
	T-Sodium	Sulfate

Instructions for Collection:	Allow sample tap to flow for a few minutes until water temperature stabilizes. Slowly fill containers to top, be sure not to over fill to prevent loss of preservative.
Required Containers/Volume:	1 L Total Chemistry – unpreserved plastic bottle 250 mL Total Metals – plastic bottle 500mL Cyanide – plastic bottle
Preservative and Handling:	Refrigerate or store on ice and do not allow to freeze, <b>recommend same day receipt at lab, must be received at lab within 24 hours of collection time</b>
Methods Technologies:	Nephelometric, Gravimetric, Flow Injection Colorimetry, Ion Chromatography, ICP, and ICPMS

Name: **Lead and Copper [ Pb ] and [ Cu ]**

Test Code: **EPA 200.8-Type 8 Metals (T-PB and T-CU)**



EPA 200.8-Type 8 Metals (Lead and Copper bottle)

Grouping:	<b>Metals</b>
Application:	Corrosion Control Assessment, Drinking Water
Analytes:	Lead (T-Pb, Total) Copper (T-Cu, Total)

Instructions for Collection:	Allow sample tap to flow for a few minutes until water temperature stabilizes. Slowly fill to line.
Required Containers/Volume:	1 L Lead and Copper – plastic bottle
Preservative and Handling:	Refrigerate or store on ice and do not allow to freeze, <b>recommend next day receipt at lab</b> , holding times, 14 days to arrive at lab for preservation with HNO <sub>3</sub> to pH <2 and 6 months to analyze
Method Technology:	Digestion technique followed by ICPMS detection

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Arsenic [ As ]**

Test Code: **EPA 200.8-T-AS , EPA 200.8-D-AS**



Grouping: **Metals**

Applications: Arsenic Rule, Drinking Water  
Surface Water

Analytes: Arsenic (T-As, Total)  
Arsenic (D-As, Dissolved)

Instructions for Collection:

Allow sample tap to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

250 mL Total Metals – plastic bottle

Preservative and Handling:


HNO<sub>3</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 6 months

Method Technology:

Digestion technique followed by ICPMS detection

Name: **Mercury [ Hg ]**

Test Codes: **EPA 245.1-T-HG , EPA 245.1-D-HG , EPA 7471B (T-HG or D-HG)**



Grouping: **Metals**

Applications: Drinking Water  
Surface Water  
Soil, solids, sludge

Analytes: Mercury (T-Hg, Total)  
Mercury (D-Hg, Dissolved)

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

Water – 250 mL Total Metals – plastic bottle  
Soils, solids, sludge – 4 oz. glass container with Teflon-lined lid

Preservative and Handling:

HNO<sub>3</sub> to pH <2, refrigerate or store on ice (4-6°C) and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days, if not preserved, 24 hours.

Method Technology:

Digestion technique followed by Cold Vapor AA detection

# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## INORGANIC AND METAL CHEMISTRY

Name: **Total-Metals and Dissolved-Metals (analytical groupings)**

Test Codes: **Total-Metals    Type 7 ,    Type 9    RCRA Types**


**Filtered-Metals    Type 3 ,    Type 4    6010-Metals (Soil)**



Type 7, RCRA, or Type 9  
(Total-Metals bottle)



Type 3 or Type 4  
(Filtered-Metals bottle)



6010-Metals  
(Soil)

### Type 3

D-Aluminum  
D-Arsenic  
D-Barium D-  
Boron  
D-Cadmium  
D-Calcium  
D-Chromium  
D-Copper  
D-Iron  
D-Lead  
D-Magnesium  
D-Manganese  
D-Mercury  
D-Nickel  
D-Potassium  
D-Selenium  
D-Silver  
D-Sodium  
D-Zinc  
Hardness

### Type 4

D-Calcium  
D-Magnesium  
D-Potassium  
D-Sodium  
Hardness

### Type 7

T-Aluminum  
T-Arsenic  
T-Barium T-  
Boron  
T-Cadmium  
T-Calcium  
T-Chromium  
T-Copper  
T-Iron  
T-Lead  
T-Manganese  
T-Magnesium  
T-Mercury  
T-Nickel  
T-Potassium  
T-Sodium  
T-Selenium  
T-Silver  
T-Zinc

### Type 9

T-Barium  
T-Cadmium  
T-Chromium  
T-Mercury  
T-Selenium

### RCRA 8

T-Arsenic T-  
Barium T-Arsenic  
T-Barium  
T-Cadmium T-  
T-Cadmium  
Chromium T-Chromium  
T-Lead T-Lead  
T-Mercury T-  
Selenium T-Mercury  
T-Silver T-Selenium  
T-Silver T-Silver

### RCRA 8+4

T-Arsenic T-  
Barium T-Arsenic  
T-Cadmium T-  
T-Cadmium  
Chromium T-Chromium  
T-Lead T-Lead  
T-Mercury T-  
Selenium T-Mercury  
T-Silver T-Selenium  
T-Silver T-Silver  
T-Copper  
T-Iron  
T-Manganese  
T-Zinc

### RCRA 8+4+6

T-Arsenic  
T-Barium  
T-Cadmium  
T-Chromium  
T-Lead  
T-Mercury  
T-Selenium  
T-Silver  
T-Copper  
T-Iron  
T-Manganese  
T-Zinc  
T-Aluminum  
T-Beryllium  
T-Cobalt  
T-Molybdenum  
T-Nickel  
T-Vanadium

**NOTE:** To test for Dissolved analytes, field filtration must be done at time of sample collection and Filtered-Metals bottle used

The constituents of the types are subject to change based upon customer request.

(Continued on Page 21)

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Total-Metals and Dissolved-Metals (analytical groupings)**

Test Codes: **Total-Metals    Type 7   ,   Type 9        RCRA Types**

**Filtered-Metals    Type 3   ,   Type 4        6010-Metals (Soil)**

**(Continued from Page 20)**

Instructions for Collection:	Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.
Required Containers/Volume:	250 mL Total Metals – plastic bottle 250 mL Filtered Metals – plastic bottle Soils, solids, sludge– 4 oz. glass container with Teflon-lined lid
Preservative and Handling:	Liquids add HNO <sub>3</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, <b>recommend next day receipt at lab</b> , holding times, 180 days to analyze, except Mercury 28 days to analyze. 24 hours for mercury if not preserved.
Methods Technologies:	Digestion techniques (except for Drinking Water) followed by analytical technique, ie: ICP (EPA 200.7, 6010), ICPMS (EPA 200.8, 6020), and Cold Vapor AA detections (EPA 245.1)

Name: **Digestion Methods for Metals (Non-drinking water)**

EPA 3010A	TCLP & Total Metal Sample Digestion
EPA 200.8	Waste water, Groundwater and Drinking Water
EPA 3050B	Solid and Hazardous Wastes

**\*NOTE: Drinking water does not require digestion unless water exceeds turbidity of 1 NTU**

#### **Analytical Methods for Metals**

EPA 200.8	Drinking water and Groundwater
EPA 200.7	Drinking water, Groundwater, and Wastewater
EPA 6010	Wastes (SW 846), Water, Soil
EPA 6020	Wastes (SW 846), Water, Soil
EPA 245.1	Mercury in Water (Prep Method 245.1)
EPA 7471B	Mercury in Soil (Prep Method 7471B)

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Ammonia [ NH<sub>3</sub> ]**

Test Code: **EPA 350.1-NH3**



Groupings: **Total-Nutrients**

Applications: Drinking Water  
Surface Water

Analyte: Ammonia

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

500 mL Nutrient Analysis – plastic bottle

Preservative and Handling:

H<sub>2</sub>SO<sub>4</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Method Technology:

Distillation technique followed by Flow Injection Colorimetry

Name: **Phosphate [ PO<sub>4</sub> ]**

Test Codes: **EPA 365.1-TPO4 , EPA 365.1-D-TP**



Groupings: **Total-Nutrients**  
**Filtered-Nutrients**

Applications: Drinking Water  
Surface Water

Analytes: Phosphate (TPO4, Total)  
Phosphate (D-TP, Dissolved)  
Phosphate (Ortho)

**NOTE:** **Prior to collection of Phosphate (Ortho), analysis must be scheduled with the laboratory at (801) 965-2400.**

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

500 mL Nutrient Analysis – plastic bottle

Preservative and Handling:

H<sub>2</sub>SO<sub>4</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days, for Phosphate (Ortho) contact laboratory

Method Technology:

Digestion technique followed by Flow Injection Colorimetry




## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Nitrate and Nitrite** [  $\text{NO}_3+\text{NO}_2$  ] , **Nitrite** [  $\text{NO}_2$  ]

Test Codes: **EPA 353.2-NO2+NO3** , **EPA 353.2-NO2**



EPA 353.2-NO2+NO3  
(Nitrate Analysis bottle)



EPA 353.2-NO2  
(Nitrite Analysis bottle)

Groupings: **Total-Nutrients**  
**Filtered-Nutrients**

Applications: Drinking Water  
Surface Water

Analytes: Nitrate/Nitrite ( $\text{NO}_2+\text{NO}_3$ )  
Nitrite ( $\text{NO}_2$ )

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

120 mL Nitrate – plastic bottle  
500 mL Total Nutrients – plastic bottle (also applicable)  
120 mL Nitrite – unpreserved plastic bottle

Preservative and Handling:

Nitrate –  $\text{H}_2\text{SO}_4$  to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

**NOTE:**

**Prior to collection of Nitrate and Nitrite, analysis must be scheduled at (801) 965-2400.**

Nitrite – No preservative, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 48 hours

Method Technology:

Flow Injection Colorimetry

Name: **Sulfate** [  $\text{SO}_4$  ]

Test Codes: **EPA 300.0-SO4C** or **EPA 375.2-SO4**



EPA 300.0-SO4C  
(Drinking Water)



EPA 375.2-SO4  
(Surface Water)

Applications: Drinking Water  
Surface Water

Analyte: Sulfate

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

120 mL Sulfate – unpreserved plastic bottle

Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Method Technology:

Ion Chromatography and Flow Injection Colorimetry

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Total-Nutrients and Dissolved-Nutrients (analytical groupings)**

Test Codes: **Total-Nutrients Type 2 , Type 3 , Type 6**

**Filtered-Nutrients Type 9**



Types 2, 3, or 6  
(Nutrient Analysis bottle)

#### Type 2

Ammonia  
Phosphate (T-PO4)

#### Type 6

Phosphate (T-PO4)  
Nitrate+Nitrite (NO2+N)

#### Type 3

Ammonia  
Phosphate (T-PO4)  
Nitrate+Nitrite (NO2+NO3)

#### Type 9

Phosphate (D-TP)  
D-Nitrate+Nitrite (NO23)  
D-Total Nitrogen



Type 9  
(Filtered-Nutrients bottle)

D = Dissolved

T = Total

**NOTE:** To test for Dissolved analytes, field filtration must be done at time of sample collection and Filtered-Nutrients bottle used.

#### Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

#### Required Containers/Volume:

500 mL Nutrient Analysis – plastic bottle  
250 mL Filtered Nutrients – plastic bottle

#### Preservative and Handling:

H<sub>2</sub>SO<sub>4</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

#### Methods Technologies:

Digestion techniques (except for Nitrate+Nitrite) followed by Flow Injection Colorimetry



## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Total-Nitrogen, Kjeldahl (TKN)**

Test Codes: **EPA 351.4**



Applications: Drinking Water  
Surface Water

Analyte: Nitrogen

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

500 mL Nutrient Analysis-plastic bottle

Preservative and Handling:

H<sub>2</sub>SO<sub>4</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Methods Technologies:

Digestion techniques followed by Flow Injection Colorimetry

Name: **Total-Nitrogen, Dissolved (DTN)**

Test Codes: **SM 4500N**



Applications: Drinking Water  
Surface Water

Analyte: Nitrogen

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

250 mL Filtered Nutrient-plastic bottle

Preservative and Handling:

H<sub>2</sub>SO<sub>4</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Methods Technologies:

Flow Injection Colorimetry


## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **pH**

Test Code: **EPA 150.1-pH**



Grouping: **Total-Chemistry**

Applications: Drinking Water  
Surface Water

Analyte: L-pH

L = Lab measured

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

1 L Total Chemistry – unpreserved plastic bottle

Preservative and Handling:


Refrigerate or store on ice and do not allow to freeze, **recommend same day receipt at lab**, total holding time 24 hours

Method Technology:

Electrometric Measurement

Name: **Conductivity**

Test Code: **EPA 120.1-COND**



Grouping: **Total-Chemistry**

Applications: Total-Metals and Filtered-Metals  
Corrosivity  
Radiologic Testing

Analyte: L-Specific Conductance

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

120 mL Conductivity – unpreserved plastic bottle  
1 L Total Chemistry – unpreserved plastic bottle (also applicable)

Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Method Technology:

Specific Conductance

Name: **Corrosivity**

Test Code: **EPA 1110 CORR**

**Note: Corrosivity requires testing of Calcium, TDS, Hardness, Alkalinity, and F-pH (field measured)**

Required Containers/Volume:

250 mL Total Metals – plastic bottle (see page 20)  
1 L Total Chemistry – plastic bottle (see page 35)

Method Technology:

Calculation: Calcium, TDS, Hardness, Alkalinity, and F-pH levels

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Color**

NOTE:

Prior to collection of Color, analysis must be scheduled at (801) 965-2400.

Test Code: **EPA 110.2-L-COLR**



Grouping: **New Drinking Water Source**

Applications: Drinking Water  
Surface Water

Analyte: Color

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume: 250 mL Color – plastic bottle

Preservative and Handling: Refrigerate or store on ice and do not allow to freeze, **recommend same day receipt at lab**, total holding time 48 hours

Method Technology: Flow Injection Colorimetry

Name: **Odor**

NOTE:

Prior to collection of Odor, analysis must be scheduled at (801) 965-2400.

Test Code: **EPA 140.1-ODOR**



Grouping: **New Drinking Water Source**

Applications: Drinking Water  
Surface Water

Analyte: Odor, Threshold Odor Number (TON)

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume: 1 L Odor – amber glass bottle

Preservative and Handling: Refrigerate or store on ice and do not allow to freeze, **recommend same day receipt at lab**, total holding time 24 hours

Method Technology: Odor Threshold (Consistent Series)

Name: **Hardness (Total Hardness as CaCO<sub>3</sub>)**

Test Code: **HARD**

Note: Hardness requires testing of Calcium and Magnesium

Required Containers/Volume: 250 mL Total Metals – plastic bottle (see page 20)

Method Technology: Calculation: Based on Calcium and Magnesium levels



## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Bromide [ Br ] , Chloride [ Cl ]**

Test Codes: **EPA 300.1-BRIC , EPA 300.0-CLIC or EPA 325.2-CL**



Grouping: **Disinfectants**

Applications: Drinking Water  
Surface Water

Analytes: Bromide  
Chloride

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

120 mL Bromide and Chloride – unpreserved plastic bottle  
120 mL Sulfate – unpreserved plastic bottle (also applicable)  
1 L Total Chemistry – unpreserved plastic bottle (also applicable)

Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Method Technology:

Ion Chromatography and Flow Injection Colorimetry

Name: **Bromate [ BrO<sub>3</sub> ] , Bromide [ Br ] , Chlorate [ ClO<sub>3</sub> ] , Chlorite [ ClO<sub>2</sub> ]**

Test Codes: **EPA 300.1-BRO3 , EPA 300.1-BRIC , EPA 300.1-CLO3 , EPA 300.1-CLO2**



Grouping: **Inorganic Disinfection By-Products**

Applications: Drinking Water  
Surface Water

Analytes: Bromate  
Bromide  
Chlorate  
Chlorite

**NOTE:** Prior to collection, analysis must be scheduled with the laboratory at (801) 965-2400

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

100 mL Bromate, Chlorate, Chlorite – plastic bottle

Preservative and Handling:

25 mg Ethylenediamine, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, holding times, 28 days to analyze, except Chlorite and Bromate require 14 days to analyze

Method Technology:

Ion Chromatography

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Perchlorate [  $\text{ClO}_4$  ]**

Test Code: **EPA 314.0-CLO4**



Application: Drinking Water

Analytes: Perchlorate  
Conductivity

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

120 mL Perchlorate – unpreserved plastic bottle  
1L Total Chemistry – unpreserved plastic bottle (also applicable)

Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze,  
**recommend next day receipt at lab**, total holding time  
28 days

Method Technology:

Ion Chromatography

Name: **Cyanide [  $\text{CN}$  ]**

Test Code: **EPA 335.4-CNCL**



Applications: Drinking Water  
Surface Water  
Groundwater

Analyte: Cyanide

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

500 mL Cyanide – plastic bottle

Preservative and Handling:

2g NaOH to pH >12, ascorbic acid in the presence of residual chlorine, refrigerate or store on ice and do not allow to freeze,  
**recommend next day receipt at lab**, total holding time 14 days

Method Technology:

Flow Injection Colorimetry (Distillation required before analysis)

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Fluoride [ F ]**

Test Code: **EPA 300.0-FLIC**



Application: Drinking Water  
Surface Water

Analyte: Fluoride

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

120 mL Fluoride – unpreserved plastic bottle

Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Method Technology:

Ion Chromatography

Name: **Sulfide [ S ]**

Test Code: **EPA 376.2-SULI**



Applications: Drinking Water  
Surface Water

Analyte: Sulfide

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

120 mL Sulfide – plastic bottle

Preservative and Handling:

3 drops Zinc Acetate and NaOH to pH >9, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 7 days

Method Technology:

Flow Injection Colorimetry




## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Total Organic Carbon (TOC)** , **UV254**

Test Codes: **SM 5310B-TOC** , **SM 5910B-UV254**



Grouping: **Total-Nutrients  
Water Treatment Technique**

Applications: Drinking Water  
Surface Water  
Groundwater

Analytes: TOC  
UV254

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

200 mL TOC – amber glass bottle  
200 mL UV254 – amber glass bottle

Preservative and Handling:

TOC – H<sub>2</sub>SO<sub>4</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

UV254 – No preservative, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 48 hours

Method Technology:

Combustion-Infrared Detection and UV Absorbance

Name: **Alkalinity**

Test Code: **SM 2320B-ALK**



Applications: Drinking Water  
Surface Water

Analytes: Alkalinity  
Carbonate Solids  
Carbonate  
Bicarbonate (BICD)  
Carbon Dioxide  
Hydroxide

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

120 mL Alkalinity – unpreserved plastic bottle  
1L Total Chemistry – unpreserved plastic bottle (also applicable)

Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 14 days

Method Technology:

pH Titration and Calculation

**INORGANIC AND METAL CHEMISTRY**

Name: **Biological Oxygen Demand (BOD)**

Test Codes: **EPA 405.1-BOD , CBOD (Carbonate buffered BOD) , SBOD (Soluble BOD) , SCBOD (Carbonations Soluble BOD)**



Applications: Drinking Water  
Surface Water

Analytes: BOD5 (5 day)  
CBOD (5 day)  
SBOD (5 day)  
SCBOD (5 day)

**Note: Please take into account that the completion of the 5 day testing duration cannot fall on a weekend or holiday.**

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

2 L BOD – unpreserved plastic bottle

Preservative and Handling:


Refrigerate or store on ice and do not allow to freeze, **recommend same day receipt at lab**, total holding time 48 hours

Method Technology:

Dissolved Oxygen Reduction over 5 days at 20 °C

Name: **Turbidity**

Test Code: **EPA 180.1-TURB**



Grouping: **Total-Chemistry  
Water Treatment Technique**

Applications: Drinking Water  
Surface Water

Analyte: Turbidity (NTU)

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

1 L Total Chemistry – unpreserved plastic bottle

Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze, **recommend same day receipt at lab**, total holding time 48 hours

Method Technology:

Nephelometric Absorbance



## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Chemical Oxygen Demand (COD)**

Test Code: **EPA 410.4-COD**



Applications: Drinking Water  
Surface Water

Analyte: COD

**Will be contracted to an outside laboratory with client's approval.**

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

500 mL Nutrient Analysis – plastic bottle

Preservative and Handling:

H<sub>2</sub>SO<sub>4</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 28 days

Method Technology:

Spectrophotometry

Name: **Chlorophyll-A**

Test Code: **SM 10200H-CH-A (modified)**



Name: **Periphyton**

Test Code: **SM 10300C (modified)**

Applications: Surface Water

Analyte: Chlorophyll-A  
Periphyton

**Note: A separate filter is required for each individual test.**

Instructions for Collection:

Chlorophyll-A: Allow source to flow for a few minutes until water temperature stabilizes. Filter 25 to 500 mL through glass fiber filter then place filter in opaque container. Also, make sure to record sample volume filtered on test request form.  
Periphyton: Follow sample collector's protocol or SM10200H (modified)

Required Containers/Volume:

Glass Fiber filter paper– store in opaque container (eg, film canister) or follow sample collector's instructions

Preservative and Handling:

Keep frozen, **recommend next day receipt at lab**, total holding time 3 weeks

Method Technology:

Homogenization followed by Spectrophotometry


## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Solids**

Test Codes: **EPA 160.1-TDS** , **EPA 160.2-TSS** , **EPA 160.4-TVS** ,  
**EPA 160.5-SS**



Applications: Drinking Water  
Surface Water

Analytes: Total Dissolved Solids (TDS) – filterable  
Total Suspended Solids (TSS) – non-filterable  
Total Volatile Solids (TVS)  
Settable Solids (SS)

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

1 L Total Chemistry – unpreserved plastic bottle

Preservative and Handling:


Refrigerate or store on ice and do not allow to freeze,  
**recommend same day receipt at lab**, holding times,  
7 days to analyze, except SS 48 hours to analyze

Method Technology:

Gravimetric detection

Name: **Silica** [ **SiO<sub>2</sub>** ]

Test Code: **EPA 370.1-SIO2**



Applications: Drinking Water  
Surface Water

Analyte: Silica (SIO2)

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

1 L Total Chemistry – unpreserved plastic bottle

Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze,  
**recommend same day receipt at lab**, total holding time 28 days

Method Technology:

Flow Injection Colorimetry


# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## INORGANIC AND METAL CHEMISTRY

Name: **Total-Chemistry (analytical groupings)**

Test Codes: **Total-Chemistry Type 2 , Type 3 , Type 9**



Types 2 or 3  
(Total-Chemistry bottle)



Type 9 (Primary Inorganics  
and Metals Chemistry)

### Type 2

pH  
Total Suspended Solids (TSS)  
Bicarbonate  
Carbon Dioxide  
Carbonate  
Chloride  
Hydroxide  
Sulfate  
Alkalinity  
Turbidity (NTU)  
L-Specific Conductance  
Total Dissolved Solids (TDS)

### Type 3

pH  
Total Suspended Solids (TSS)  
D-Calcium  
D-Magnesium  
D-Potassium  
D-Sodium  
Bicarbonate  
Carbon Dioxide  
Carbonate  
Chloride  
Hydroxide  
Sulfate  
Alkalinity  
Hardness  
Turbidity (NTU)  
L-Specific Conductance  
Total Dissolved Solids (TDS)  
Carbonate Solids

### Type 9

Cyanide  
Fluoride  
Turbidity (NTU)  
Total Dissolved Solids (TDS)  
T-Arsenic  
T-Barium  
T-Beryllium  
T-Cadmium  
T-Chromium  
T-Copper  
T-Lead  
T-Mercury  
T-Nickel  
T-Selenium  
T-Antimony  
T-Thallium  
T-Sodium  
Sulfate

D = Dissolved

T = Total

L = Lab measured

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

1 L Total Chemistry – unpreserved plastic bottle

Type 9 also under Annual Inorganics and Metals (see page 19)

1 L Total Chemistry – unpreserved plastic bottle

250 mL Total Metals – plastic bottle

500mL Cyanide – plastic bottle

Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, total holding time 48 hours

Methods Technologies:

Electrometric, Nephelometric, Gravimetric, Titration, Flow Injection Colorimetry, Ion Chromatography, ICP, and ICPMS


## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### INORGANIC AND METAL CHEMISTRY

Name: **Chromium-VI [ Cr<sup>6+</sup> ]**

Test Code: **EPA 218.7**



Analyte: Chromium-VI (Hexavalent)

**NOTE:** Prior to collection of Chromium-VI samples, analysis must be scheduled at (801) 965-2400.

Instructions for Collection:	Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.
Required Containers/Volume:	120 mL Chromium-VI – plastic bottle
Preservative and Handling:	1 mL of (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> /NH <sub>4</sub> OH to 100 mL. Refrigerate or store on ice and do not allow to freeze, <b>recommend same day receipt at lab</b> , total holding time 14 days
Method Technology:	Ion Chromatography

Name: **Uranium [ U ]**

Test Codes: **EPA 200.8-UUMS (Total)** , **EPA 200.8-UFMS (Dissolved)**



Applications: Drinking Water  
Surface Water

Analytes: T-Thallium  
T-Uranium (Unfiltered)  
D-Uranium (Filtered)

Instructions for Collection:	Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill containers to top, be sure not to over fill to prevent loss of preservative.
Required Containers/Volume:	250 mL Total Metals – plastic bottle 2 L RadioChem – plastic bottle (also applicable for Uranium)
Preservative and Handling:	HNO <sub>3</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, <b>recommend next day receipt at lab</b> , total holding time 180 days
Method Technology:	Digestion techniques (except for Drinking Water) followed by ICPMS detection


# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## INORGANIC AND METAL CHEMISTRY

Name: **RCRA Toxic Characteristic Leaching Procedure (TCLP)**

Test Code: **TCLP- Metals**



Grouping: **Inorganic Testing (RCRA)**

**NOTE:** Prior to collection of 1311-TCLP (water and soil) or 1332-TCLP (oil), analysis must be scheduled at (801) 965-2400.

Analytes:

### RCRA 8

D004 T-Arsenic  
D005 T-Barium  
D006 T-Cadmium  
D007 T-Chromium  
D008 T-Lead  
D009 T-Mercury  
D010 T-Selenium  
D011 T-Silver

### RCRA 8+4

D004 T-Arsenic  
D005 T-Barium  
D006 T-Cadmium  
D007 T-Chromium  
D008 T-Lead  
D009 T-Mercury  
D010 T-Selenium  
D011 T-Silver  
T-Copper  
T-Iron  
T-Manganese  
T-Zinc

### RCRA 8+4+6

D004 T-Arsenic  
D005 T-Barium  
D006 T-Cadmium  
D007 T-Chromium  
D008 T-Lead  
D009 T-Mercury  
D010 T-Selenium  
D011 T-Silver  
T-Copper  
T-Iron  
T-Manganese  
T-Zinc  
T-Aluminum  
T-Beryllium  
T-Cobalt  
T-Molybdenum  
T-Nickel  
T-Vanadium

\*Other metals may be analyzed, but must be specified on test request form. For complete list of metals see indexes or contact us 801-965-2400.

Instructions for Collection:

Slowly fill to top of container.

Required Containers/Volume:

Water - 4 L amber glass bottle

Soil - 4 oz. amber with Teflon-lined lid glass container

Preservative and Handling:

Refrigerate or store on ice 4-6°C and do not allow to freeze, **recommend next day receipt at lab**, holding times, Metals - 7 days to TCLP and 180 days to analyze except Mercury - 7 days to TCLP and 28 days to analyze

Method Technology:

Method 1311 Leaching procedure and extraction technique followed by Method 6010 ICP, Method 6020 ICPMS, and Cold Vapor AA detections

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### ENVIRONMENTAL MICROBIOLOGY

Name: **Total Coliform** , **Fecal Coliform** and ***E. coli***

Test Codes: **SM 9223B-Colilert** , **SM 9222B/9221E-mENDO+EC** ,

**SM9222B/9222D-mENDO+mFC**



SM 9223B-Colilert  
SM 9222B/9221E-  
mENDO+EC  
(Drinking Water)



SM 9222B/9222D-  
mENDO+mFC  
(Surface Water)

Applications: Drinking Water  
Surface Water

Analytes: Total Coliforms  
Fecal Coliforms  
*E. coli* (Colilert test only)

**NOTE: For Surface Water: Sample delivery should be within 2 hours and temperature should not exceed 10°C.**

Instructions for Collection:	Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to between the line and the top of the bottle, be sure not to over fill to prevent loss of preservative.
Required Containers/Volume:	120 mL Water Bacteriology – sterile plastic bottle
Preservative and Handling:	10 mg Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> , refrigerate or store on ice and do not allow to freeze, <b>recommend same day receipt at lab</b> , total holding time 30 hours, except for Surface Water samples 8 hours
Method Technology:	Chromofluorogenic, Membrane Filtration, and Fermentation

---

Name: **Heterotrophic Plate Count (HPC)**

Test Code: **SM 9215B-HPC**



SM 9215B-HPC  
(Water Bacteriology bottle)

Applications: Drinking Water  
Pool/Spa Water  
Reagent Water (Deionized, RO, Distilled, etc.)  
Surface Water

Analyte: Heterotrophic Plate Count (total bacteria count)

Instructions for Collection:	Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.
Required Containers/Volume:	120 mL Water Bacteriology – sterile plastic bottle
Preservative and Handling:	10 mg Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> , refrigerate or store on ice and do not allow to freeze, <b>recommend same day receipt at lab</b> , total holding time 30 hours, except for Surface Water samples 8 hours
Method Technology:	Pour Plate Agar



## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory


### ENVIRONMENTAL MICROBIOLOGY

Name: ***Cryptosporidium* and *Giardia***

Test Code: **EPA 1623-Protozoa**



EPA 1623-Protozoa Filter  
(Pall *Envirochek HV*)



EPA 1623-Protozoa Matrix  
Spike  
(20 L Carboy)

Grouping: **Water Treatment Technique**

Application: Surface Water Treatment Rule (LT2)

Analytes: *Cryptosporidium* (oocysts)  
*Giardia* (cysts)

**Will be contracted to an outside laboratory with client's approval.**

#### Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. For bulk water samples, slowly fill to top of container. For field filtration, pressurized flow is necessary to pass the required 10 L (2.6 gal) volume through the capsule filter (a pump may be needed for stream flow analysis or where pressurized flow is not available). Also, make sure to record sample volume filtered on test request form.

#### Required Containers/Volume:

1  $\mu$ m pore size Pall *Envirochek HV* – capsule filter  
20 L Carboy – plastic bottle (6-1 gal plastic jugs also applicable)

#### Preservative and Handling:

Refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, holding times, 96 hours to elution, 72 hours to staining, and 7 days to complete microscopy analysis

For regulatory analysis, samples must arrive at laboratory at less than 10 °C or arrive the same day as collected

#### Method Technologies:

Filtration, Elution, Centrifugation, Immunomagnetic Separation, Fluorescent-conjugated Antibody Staining, and Advanced Microscopy with Differential Interference Contrast

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**ENVIRONMENTAL MICROBIOLOGY**

Name: ***Legionella***

Test Code: **SM 9260J-Legionella**



Applications: Drinking Water (Hot Water Tanks)  
Air Handling (Swamp Coolers, Evaporators, etc.)

Analyte: *Legionella*

**NOTE:** Prior to collection of SM 9260J-Legionella samples, analysis must be scheduled at (801) 965 -2400. Batch analysis performed on third Thursday of each month).

Instructions for Collection: Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume: 2 L Legionella – unpreserved plastic bottle

Preservative and Handling: **Same day receipt at lab**, lab testing is started within 48 hours of receiving sample

Method Technology: Filtration, Elution, MWY Agar, and Latex-Agglutination

**\*\*The State of Utah Public Health Laboratory is a member of the CDC's elite program for Legionella testing.**

Name: ***Enterococcus***

Test Code: **Type 9-ECOCCL**

**Will be contracted to an outside laboratory with client's approval.**

**Note:** *Enterococcus* tested with Total/Fecal Coliforms for Surface Water analysis

Required Containers/Volume: 120 mL Water Bacteriology – sterile plastic bottle (see page 38)

Method Technology: Enzyme Substrate



## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### RADIOLOGIC CHEMISTRY

Name: **Gross-Alpha , Gross-Beta , Uranium [ U ]**

Test Codes: **EPA 900.0-ALPG , EPA 900.0-BETG , EPA 200.8-UUMS\* ,  
EPA 200.8-UFMS\***



RadioChem & Conductivity  
(Drinking and Surface Water)

Grouping: **Radiologic Testing**

Applications: Drinking Water  
Surface Water

Analytes: Gross-Alpha  
Gross-Beta  
Uranium (UUMS-Unfiltered)  
Uranium (UFMS-Filtered)  
Conductivity

**NOTE: If Gross-alpha level is above 5 pCi/l, Radium-226 (page 42) must also be tested.**

**Uranium also listed under Inorganic Testing (see page 36).**

#### Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative. Do not pour from RadioChem bottle to Conductivity bottle, must fill Conductivity separately.

#### Required Containers/Volume:

##### Drinking Water

2 L RadioChem – plastic bottle

\*250 mL Total Metals – plastic bottle (also applicable for Uranium)

120 mL Conductivity – plastic bottle

##### Surface Water

2 L RadioChem – plastic bottles

\*250 mL Total Metals – plastic bottle (also applicable for Uranium)

120 mL Conductivity – plastic bottle

#### Preservative and Handling:

HNO<sub>3</sub> to pH <2 (No preservation for conductivity bottle), **recommend next day receipt at lab**, total holding time 180 days

#### Methods Technologies:

Evaporation, Co-Precipitation, and Gas-Flow Proportional Counter  
Uranium detected by ICPMS

## **RADIOLOGIC CHEMISTRY**

Name: **Radium-226 , Radium-228**

Test Codes: **EPA 903.1-226R , EPA 904.0-228R**



Grouping: **Radiologic Testing**

Applications: Drinking Water  
Surface Water

Analytes: Radium-226  
Radium-228

Instructions for Collection:	Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.
Required Containers/Volume:	<u>Drinking Water</u> 2 L RadioChem – plastic bottle  <u>Surface Water</u> 2 L RadioChem – plastic bottles
Preservative and Handling:	HNO <sub>3</sub> to pH <2, refrigerate or store on ice and do not allow to freeze, <b>recommend next day receipt at lab</b> , total holding time 180 days
Methods Technologies:	Evaporation, Co-Precipitation, Emanation, and Gas-Flow Proportional Counter

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### RADIOLOGIC CHEMISTRY

Name: **Radon [ Rn ] , Gross Gamma**

Test Codes: **EPA 913.0-RN , EPA 901.1-GAMG**




Grouping: **Radiologic Testing**

Applications: Drinking Water  
Surface Water

Analytes: Radon  
Gamma Emission  
Conductivity



**NOTE:** Radon and Gross Gamma testing will be contracted to an outside laboratory with client's approval.



#### Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative. For Radon, check for air bubbles by inverting and fill the remaining portion if bubbles found. There must be no headspace. Also, do not pour from RadioChem bottle to Conductivity bottle, must fill Conductivity separately.

#### Required Containers/Volume:

##### Radon

3-40 mL Radon – unpreserved glass vials

##### Gross Gamma (Drinking Water)

2 L RadioChem – plastic bottle

120 mL Conductivity – plastic bottle

##### Gross Gamma (Surface Water)

2-2 L RadioChem – plastic bottles

120 mL Conductivity – plastic bottle

#### Preservative and Handling:

Radon – No preservative, use insulated packaging, **mandatory next day receipt at lab**, total holding time 72 hours

Gross Gamma – HNO<sub>3</sub> to pH <2, **recommend next day receipt at lab**, total holding time 180 days

#### Methods Technologies:

LS Counter and Gamma Spectrometer

## ***Chemical and Environmental Laboratory – Client Services Manual***

State of Utah Public Health Laboratory

### **SPECIALTY ANALYSES**

Name: **Geosmin and MIB**

Test Code: **EPA 525.2-ODOR**



Applications: Drinking Water

Analytes: Geosmin  
2-Methylisoborneol (MIB)

**NOTE:** Prior to collection of EPA 525.2-ODOR samples, analysis must be scheduled at (801) 965-2400, due to complex nature of testing procedures (batch analysis is preferred).

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container.

Required Containers/Volume:

1 L Method EPA 525.2 Odor – amber glass bottle

Preservative and Handling:

**Same day receipt at lab**, no holding time

Method Technology:

Liquid-solid extraction followed by GCMS detection

Name: **Lead (air samples)**

Test Code: **EPA 200.8-T-Pb**



Application: Air Quality

Analyte: Lead (air sample particulate analysis)

Instructions for Collection:

Air monitoring technique to capture particulates using filter paper.

Required Containers/Volume:

Preservative and Handling:

Store filtered samples in plastic bag

Method Technology:

Filter digestion followed by ICPMS detection

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### SPECIALTY ANALYSES

#### Unregulated Contaminant Monitoring Rule3 (UCMR3) Project (2013-2015)

According to the 1996 Safe Drinking Water Act (SDWA) a new list of 30 unregulated contaminants must be monitored every 5 years. The Chemical and Environmental Laboratory is certified by the EPA to test UCMR 1,2, and now 3 (List 1) from January 2013 through December 2015.

- Attn: Open tap and run 5 min before taking sample.
- Refrigerate all samples
- If shipping, add appropriate coolant
- Additional 537 & 522 bottles for Matrix spike/Matrix spike duplicate

- One extra bottle for matrix spike for EPTDS only



List 1/EPTDS-16 Bottles

Grouping: List 1  
Tests (6 Total):

EPA 537	EPTDS
EPA 524.3	EPTDS
EPA 522	EPTDS
EPA 200.8	EPTDS & DSMRT
EPA 218.7	EPTDS & DSMRT
EPA 300.1	EPTDS & DSMRT

Name: **UCMR3 Perfluorinated Compounds**

Test Codes: **EPA 537—EPTDS**



Sample Bottles Trip Blanks

Instructions for Collection:

- Fill 2 Sample Bottles
- Transfer sealed Trip Blank to the empty Blank Bottle

Containers/Volume:

Sample: Slowly fill bottle to top of container, be sure not to over fill to prevent loss of preservative. Agitate to dissolve preservative. Field blank: Pour all contents of Trip Blank into Blank bottle at test site. Discard Trip Blank bottle.

4/250 mL Plastic; 2 Sample, 1 Blank, 1 Trip Blank

Name: **UCMR3 Volatiles**

Test Codes: **EPA 524.3—EPTDS**



Sample Bottles

Trip Blanks

Instructions for Collection:

- Do not open Trip Blanks
- Fill 4 Sample Vials with no headspace

Containers/Volume:

Slowly fill 4 Sample vials to top of container, be sure not to over fill. Remove air bubbles by inverting, and fill the remaining portion if bubbles found. There must be no headspace.

6/40 mL clear VOAs; 4 Sample; 2 Trip Blanks

## Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

### SPECIALTY ANALYSES

Name: **UCMR3 1,4 Dioxane**

Test Codes: **EPA 522—EPTDS**



Instructions for Collection:

- **Fill 2 Sample Bottles**
- **Agitate to dissolve**
- **Add 4 mL vial, cap and mix**

Slowly fill sample bottle to neck, be sure not to over fill to prevent loss of preservative. Agitate to dissolve preservative. Once dissolved, add contents of 4 mL vial, cap, and mix.

Containers/Volume:

2/500 mL AG & 2/4 mL vial

Name: **UCMR3 Metals**

Test Codes: **EPA 200.8—EPTDS & DSMRT**



Instructions for Collection:

- **Uncap sealed Trip Blank**
- **Fill the Sample Bottle**
- **Recap the Trip Blank**
- **List Location on Trip Blank**

Remove cap from blank exposure and place it facedown on a clean flat surface. Wait 5 seconds. Slowly fill sample bottle to top of container, be sure not to over fill to prevent loss of preservative. Agitate to dissolve preservative. Reseal blank bottle.

Containers/Volume:

2/250 mL Plastic

Name: **UCMR3 Hexavalent Chromium**

Test Codes: **EPA 218.7—EPTDS & DSMRT**



Instructions for Collection:

- **Fill 1 Sample Bottle**

Slowly fill sample bottle to top of container, be sure not to over fill to prevent loss of preservative. Agitate to dissolve preservative. Order new bottle if preservative expiration date has past.

Containers/Volume:

100 mL Plastic

Name: **UCMR3 Chlorate**

Test Codes: **EPA 300.1—EPTDS & DSMRT**



Instructions for Collection:

- **Fill 1 Sample Bottle**

Slowly fill bottle to top of container, be sure not to over fill to prevent loss of preservative. Agitate to dissolve preservative. (If using chlorine dioxide, sparge with an inert gas for 5 min before preservation).

Containers/Volume:

100 mL Plastic



**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Clean Water Act (CWA)**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Alkalinity	<b>SM 2320B-ALK</b>	Inorganic Chemistry	<b>30,34</b>
Aluminum	<b>EPA 200.8-AL</b>	Inorganic Chemistry-Metals	<b>20</b>
Ammonia	<b>EPA 350.1-NH3</b>	Inorganic Chemistry-Nutrients	<b>22,24</b>
Antimony	<b>EPA 200.8-SB</b>	Inorganic Chemistry-Metals	<b>20</b>
Arsenic	<b>EPA 200.8-AS</b>	Inorganic Chemistry-Metals	<b>19,20</b>
Barium	<b>EPA 200.8-BA</b>	Inorganic Chemistry-Metals	<b>20</b>
Beryllium	<b>EPA 200.8-BE</b>	Inorganic Chemistry-Metals	<b>20</b>
BOD	<b>EPA 405.1-BOD</b>	Inorganic Chemistry	<b>32</b>
Boron	<b>EPA 200.7-B</b>	Inorganic Chemistry-Metals	<b>20</b>
Cadmium	<b>EPA 200.8-CD</b>	Inorganic Chemistry-Metals	<b>20</b>
Calcium	<b>EPA 200.7-CA</b>	Inorganic Chemistry-Metals	<b>20</b>
CBOD	<b>CBOD</b>	Inorganic Chemistry	<b>32</b>
Chloride	<b>EPA 325.2-CL</b>	Inorganic Chemistry	<b>28,35</b>
Chlorophyll-A	<b>SM 10200H-CH-A</b>	Inorganic Chemistry	<b>33</b>
Chromium	<b>EPA 200.8-CR</b>	Inorganic Chemistry-Metals	<b>20</b>
Chromium-IV	<b>EPA 218.7</b>	Inorganic Chemistry	<b>36,46</b>
Cobalt	<b>EPA 200.8-CO</b>	Inorganic Chemistry-Metals	<b>20</b>
COD	<b>EPA 410.4-COD</b>	Inorganic Chemistry	<b>32</b>
Conductivity	<b>EPA 120.1-COND</b>	Inorganic Chemistry	<b>26,41</b>
Copper	<b>EPA 200.8-CU</b>	Inorganic Chemistry-Metals	<b>20</b>
Cyanide	<b>EPA -CNCL</b>	Inorganic Chemistry	<b>29,35</b>
<i>E. coli</i>	<b>SM 9223B-Colilert</b>	Environmental Microbiology	<b>38</b>
<i>Enterococcus</i>	<b>Enterolert</b>	Environmental Microbiology	<b>40</b>
Fecal Coliform	<b>SM 9222D-mFC</b>	Environmental Microbiology	<b>38</b>
Fluoride	<b>EPA 300.0-FLIC</b>	Inorganic Chemistry	<b>30,35</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Clean Water Act (CWA)**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Gross Alpha	<b>EPA 900.0-ALPG</b>	Radiologic Chemistry	<b>41</b>
Gross Beta	<b>EPA 900.0-BETG</b>	Radiologic Chemistry	<b>41</b>
Gross Gamma	<b>EPA 901.1-GAMG</b>	Radiologic Chemistry	<b>43</b>
Hardness	<b>HARD</b>	Inorganic Chemistry	<b>27,35</b>
Heterotrophic Plate Count	<b>SM 9215B-HPC</b>	Environmental Microbiology	<b>38</b>
Iron	<b>EPA 200.7-FE</b>	Inorganic Chemistry-Metals	<b>20</b>
Lead	<b>EPA 200.8-PB</b>	Inorganic Chemistry-Metals	<b>20</b>
Magnesium	<b>EPA 200.7-MG</b>	Inorganic Chemistry-Metals	<b>20</b>
Manganese	<b>EPA 200.8-MN</b>	Inorganic Chemistry-Metals	<b>20</b>
Mercury	<b>EPA 245.1-HG</b>	Inorganic Chemistry-Metals	<b>19,20</b>
Mercury	<b>EPA 200.8-HG</b>	Inorganic Chemistry-Metals	<b>19,20</b>
Molybdenum	<b>EPA 200.8-MO</b>	Inorganic Chemistry-Metals	<b>20</b>
Nickel	<b>EPA 200.8-NI</b>	Inorganic Chemistry-Metals	<b>20</b>
Nitrate+Nitrite	<b>EPA 353.2-NO2+NO3</b>	Inorganic Chemistry-Nutrients	<b>23,24</b>
Nitrite (only)	<b>EPA 353.2-NO2</b>	Inorganic Chemistry-Nutrients	<b>23,24</b>
Oil and Grease	<b>EPA 1664-O/G</b>	Organic Chemistry	<b>14</b>
Phosphate	<b>EPA 365.1-TPO4</b>	Inorganic Chemistry-Nutrients	<b>22,24</b>
pH	<b>EPA 150.1-PH</b>	Inorganic Chemistry	<b>25,34</b>
Potassium	<b>EPA 200.7-K</b>	Inorganic Chemistry-Metals	<b>20</b>
Radium-226	<b>EPA 903.1-226R</b>	Radiologic Chemistry	<b>42</b>
Radium-228	<b>EPA 904.0-226R</b>	Radiologic Chemistry	<b>42</b>
Radon	<b>EPA 913.0-RN</b>	Radiologic Chemistry	<b>43</b>
Selenium	<b>EPA 200.8-SE</b>	Inorganic Chemistry-Metals	<b>20</b>
Silica	<b>EPA 370.1-SIO2</b>	Inorganic Chemistry	<b>34</b>
Silver	<b>EPA 200.8-AG</b>	Inorganic Chemistry-Metals	<b>20</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Clean Water Act (CWA)  
Clean Water Act (CWA) – EPA 608-PCB/OcPest**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Sodium	<b>EPA 200.7-NA</b>	Inorganic Chemistry-Metals	<b>20</b>
Sulfate	<b>EPA 375.2-SO4</b>	Inorganic Chemistry	<b>23,35</b>
Sulfide	<b>EPA 376.2-SULI</b>	Inorganic Chemistry	<b>30</b>
Thallium	<b>EPA 200.8-TL</b>	Inorganic Chemistry-Metals	<b>35</b>
TOC	<b>SM 5310B-TOC</b>	Inorganic Chemistry	<b>30</b>
Total Coliform	<b>SM 9223B-Colilert</b>	Environmental Microbiology	<b>38</b>
Total Coliform	<b>SM 9222B-mENDO</b>	Environmental Microbiology	<b>38</b>
Total Dissolved Solids (TDS)	<b>EPA 160.1-TDS</b>	Inorganic Chemistry	<b>34,35</b>
Total Suspended Solids (TSS)	<b>EPA 160.2-TSS</b>	Inorganic Chemistry	<b>34,35</b>
Turbidity	<b>EPA EPA 180.1-TURB</b>	Inorganic Chemistry	<b>32,35</b>
Uranium	<b>EPA 200.8-UUMS</b>	Inorganic Chemistry-Metals	<b>36,41</b>
UV254	<b>SM 5910B-UV254</b>	Inorganic Chemistry	<b>31</b>
Vanadium	<b>EPA 200.8-V</b>	Inorganic Chemistry-Metals	<b>20</b>
Zinc	<b>EPA 200.8-ZN</b>	Inorganic Chemistry-Metals	<b>20</b>
4,4'-DDD	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
4,4'-DDE	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
4,4'-DDT	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Aldrin	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
alpha-BHC	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
beta-BHC	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Chlordane	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
delta-BHC	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Dieldrin	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Endosulfan I	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>

# **Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

## **PROGRAM SPECIFIC INDEXES**

### **Clean Water Act (CWA) – EPA 608-PCB/OcPest Clean Water Act (CWA) – EPA 624-VOC**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Endosulfan II	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Endosulfan sulfate	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Endrin	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Endrin aldehyde	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
gamma-BHC	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Heptachlor	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Heptachlor epoxide	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Methoxychlor	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
PCB-1016	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
PCB-1221	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
PCB-1232	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
PCB-1242	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
PCB-1248	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
PCB-1254	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
PCB-1260	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
Toxaphene	<b>EPA 608-PCB/OcPest</b>	Organic Chemistry-PCB/OcPest	<b>14</b>
1,1-Dichloroethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,1-Dichloroethene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,1-Dichloropropene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,1,1-Trichloroethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,1,1,2-Tetrachloroethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,1,2-Trichloroethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,1,2-Trichloro-1,2,2-trifluoroethane (Freon 113)	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,1,2,2-Tetrachloroethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>

# **Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

## **PROGRAM SPECIFIC INDEXES**

### **Clean Water Act (CWA) – EPA 624-VOC**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
1,2-Dibromo-3-chloropropane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,2-Dichlorobenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,2-Dichloroethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,2-Dichloropropane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,2-Dichlorotoluene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,2,3-Trichlorobenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,2,3-Trichloropropane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,2,4-Trichlorobenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,2,4-Trimethylbenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,3-Dichlorobenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,3-Dichloropropane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,3,5-Trimethylbenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,4-Dichlorobenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,4-Dichlorotoluene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,4-Isopropyltoluene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
2,2-Dichloropropane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Benzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Bromobenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Bromochloromethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Bromodichloromethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Bromoform	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Bromomethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Carbon tetrachloride	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Chlorobenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Chlorodibromomethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>

# ***Chemical and Environmental Laboratory – Client Services Manual***

State of Utah Public Health Laboratory

## **PROGRAM SPECIFIC INDEXES**

### **Clean Water Act (CWA) – EPA 624-VOC**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Chloroethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Chloroform	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Chloromethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
cis-1,2-Dichloroethene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
cis-1,3-Dichloropropene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Dibromomethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Dichlorodifluoromethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Ethylbenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Ethylene dibromide	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Hexachlorobutadiene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Isopropylbenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Methylene chloride	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
MTBE	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Napthalene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
n-Butylbenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
n-Propylbenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
sec-Butylbenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Styrene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
tert-Butylbenzene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Tetrachloroethene (PCE)	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Toluene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
trans-1,2-Dichloroethene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
trans-1,3-Dichloropropene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Trichloroethene (TCE)	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Trichlorofluoromethane	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>



# **Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

## **PROGRAM SPECIFIC INDEXES**

### **Clean Water Act (CWA) – EPA 624-VOC Clean Water Act (CWA) – EPA 625-SVOC**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Vinyl chloride	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
Xylene	<b>EPA 624-VOC</b>	Organic Chemistry-VOC	<b>9</b>
1,2-Dichlorobenzene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
1,2,4-Trichlorobenzene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
1,3-Dichlorobenzene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
1,4-Dichlorobenzene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2-Chloronaphthalene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2-Chlorophenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2-Methyl naphthalene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2-Methyl phenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2-Methyl-4,6-dinitrophenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2-Nitroaniline	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2-Nitrophenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2,4-Dichlorophenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2,4-Dimethylphenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2,4-Dinitrophenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2,4-Dinitrotoluene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2,4,6-Trichlorophenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
2,6-Dinitrotoluene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
3-Methyl phenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
3-Nitroaniline	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
3,3'-Dichlorobenzidine	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
4-Bromophenyl phenyl ether	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
4-Chloroaniline	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>

# **Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

## **PROGRAM SPECIFIC INDEXES**

### **Clean Water Act (CWA) – EPA 625-SVOC**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
4-Chlorophenyl phenyl ether	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
4-Chloro-3-methyl phenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
4-Methyl phenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
4-Nitroaniline	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
4-Nitrophenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Acenaphthene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Acenaphthylene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Aniline	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Anthracene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Benzdine	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Benzo (a) anthracene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Benzo (a) pyrene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Benzo (b) fluoranthene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Benzo (g,h,i) perylene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Benzo (k) fluoranthene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Benzylbutylphthalate	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Benzyl alcohol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Benzoic acid	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
bis (2-chloroethyl) ether	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
bis (2-chloroethoxy) methane	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
bis (2-chloroisopropyl) ether	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
bis (2-ethylhexyl) phthalate	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Chrysene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Dibenzo (a,h) anthracene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Dibenzofuran	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES**

**Clean Water Act (CWA) – EPA 625-SVOC**

Analyte	Test Code	UPHL Unit	Page Number
Diethyl phthalate	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Dimethyl phthalate	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Di-n-butyl phthalate	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Di-n-octyl phthalate	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Fluoranthene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Fluorene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Hexachlorobenzene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Hexachlorobutadiene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Hexachlorocyclopentadiene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Hexachloroethane	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Indeno (1,2,3-cd) pyrene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Isophorone	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
n-Nitrosodimethylamine	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
n-Nitrosodiphenylamine	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
n-Nitrosodipropylamine	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Naphthalene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Nitrobenzene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Pentachlorophenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Phenanthrene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Phenol	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>
Pyrene	<b>EPA 625-SVOC</b>	Organic Chemistry-SVOC	<b>10</b>

# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## PROGRAM SPECIFIC INDEXES

### Resource Conservation and Recovery Act (RCRA) – EPA 8081-PCB/OcPest

Analyte	Test Code	EPA HW No.	UPHL Unit	Page Number
4,4'-DDD	EPA 8081-OcPest	U060	Organic Chemistry-PCB/OcPest	14
4,4'-DDE	EPA 8081-OcPest		Organic Chemistry-PCB/OcPest	14
4,4'-DDT	EPA 8081-OcPest	U061	Organic Chemistry-PCB/OcPest	14
Aldrin	EPA 8081-OcPest	P005	Organic Chemistry-PCB/OcPest	14
alpha-BHC	EPA 8081-OcPest	U129	Organic Chemistry-PCB/OcPest	14
beta-BHC	EPA 8081-OcPest		Organic Chemistry-PCB/OcPest	14
Chlordane	EPA 8081-OcPest	D020	Organic Chemistry-PCB/OcPest	14,16
delta-BHC	EPA 8081-OcPest		Organic Chemistry-PCB/OcPest	14
Dieldrin	EPA 8081-OcPest	P037	Organic Chemistry-PCB/OcPest	14
Endosulfan I	EPA 8081-OcPest	P050	Organic Chemistry-PCB/OcPest	14
Endosulfan II	EPA 8081-OcPest	P050	Organic Chemistry-PCB/OcPest	14
Endosulfan sulfate	EPA 8081-OcPest	P050	Organic Chemistry-PCB/OcPest	14
Endrin	EPA 8081-OcPest	D012	Organic Chemistry-PCB/OcPest	14,16
Endrin aldehyde	EPA 8081-OcPest	D012	Organic Chemistry-PCB/OcPest	14
gamma-BHC (Lindane)	EPA 8081-OcPest	D013	Organic Chemistry-PCB/OcPest	14,16
Heptachlor	EPA 8081-OcPest	D031	Organic Chemistry-PCB/OcPest	14,16
Heptachlor epoxide	EPA 8081-OcPest	D031	Organic Chemistry-PCB/OcPest	14,16
Methoxychlor	EPA 8081-OcPest	D014	Organic Chemistry-PCB/OcPest	14,16
PCB-1016	EPA 8082-PCB		Organic Chemistry-PCB/OcPest	14
PCB-1221	EPA 8082-PCB		Organic Chemistry-PCB/OcPest	14
PCB-1232	EPA 8082-PCB		Organic Chemistry-PCB/OcPest	14
PCB-1242	EPA 8082-PCB		Organic Chemistry-PCB/OcPest	14
PCB-1248	EPA 8082-PCB		Organic Chemistry-PCB/OcPest	14
PCB-1254	EPA 8082-PCB		Organic Chemistry-PCB/OcPest	14
PCB-1260	EPA 8082-PCB		Organic Chemistry-PCB/OcPest	14

# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## PROGRAM SPECIFIC INDEXES

### Resource Conservation and Recovery Act (RCRA) – EPA 8081-PCB/OcPest Resource Conservation and Recovery Act (RCRA) – EPA 8151-Herbicides

Analyte	Test Code	EPA HW No.	UPHL Unit	Page Number
Toxaphene	EPA 8081-OcPest	D015	Organic Chemistry- PCB/OcPest	14,16
*Add-ons	EPA 8081-OcPest		Organic Chemistry- PCB/OcPest	14
2,4-D	EPA 8151-Herbicides	D016	Organic Chemistry-Herbicides	11,16
2,4,5-TP (Silvex)	EPA 8151-Herbicides	D017	Organic Chemistry-Herbicides	11,16
Dalapon	EPA 8151-Herbicides		Organic Chemistry-Herbicides	11
Dinoseb	EPA 8151-Herbicides	P020	Organic Chemistry-Herbicides	11
Pentachlorophenol	EPA 8151-Herbicides	D037	Organic Chemistry-Herbicides	11,16
Picloram	EPA 8151-Herbicides		Organic Chemistry-Herbicides	11
*Add-ons	EPA 8151-Herbicides		Organic Chemistry-Herbicides	11
1,1-Dichloroethane	EPA 8260-VOC		Organic Chemistry-VOC	9
1,1-Dichloroethene	EPA 8260-VOC	D029	Organic Chemistry-VOC	9,16
1,1-Dichloropropene	EPA 8260-VOC		Organic Chemistry-VOC	9
1,1,1-Trichloroethane	EPA 8260-VOC		Organic Chemistry-VOC	9
1,1,1,2-Tetrachloroethane	EPA 8260-VOC	U208	Organic Chemistry-VOC	9
1,1,2-Trichloroethane	EPA 8260-VOC	U227	Organic Chemistry-VOC	9
1,1,2-Trichloro-1,2,2-trifluoroethane (Freon 113)	EPA 8260-VOC		Organic Chemistry-VOC	9
1,1,2,2-Tetrachloroethane	EPA 8260-VOC	U209	Organic Chemistry-VOC	9
1,2-Dibromo-3-chloropropane	EPA 8260-VOC	U066	Organic Chemistry-VOC	9
1,2-Dichlorobenzene	EPA 8260-VOC	U070	Organic Chemistry-VOC	9
1,2-Dichloroethane	EPA 8260-VOC	D028	Organic Chemistry-VOC	9,16
1,2-Dichloropropane	EPA 8260-VOC	U083	Organic Chemistry-VOC	9
1,2-Dichlorotoluene	EPA 8260-VOC		Organic Chemistry-VOC	9
1,2,3-Trichlorobenzene	EPA 8260-VOC		Organic Chemistry-VOC	9

\*Tentative analysis by special request, record specific analyte on test request form and schedule at 801-965-2400

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Resource Conservation and Recovery Act (RCRA) – EPA 8260-VOC**

Analyte	Test Code	EPA HW No.	UPHL Unit	Page Number
1,2,3-Trichloropropane	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
1,2,4-Trichlorobenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
1,2,4-Trimethylbenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
1,3-Dichlorobenzene	<b>EPA 8260-VOC</b>	U071	Organic Chemistry-VOC	<b>9</b>
1,3-Dichloropropane	<b>EPA 8260-VOC</b>	U084	Organic Chemistry-VOC	<b>9</b>
1,3,5-Trimethylbenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
1,4-Dichlorobenzene	<b>EPA 8260-VOC</b>	D027	Organic Chemistry-VOC	<b>9,16</b>
1,4-Dichlorotoluene	<b>EPA 8260-VOC</b>	U072	Organic Chemistry-VOC	<b>9</b>
1,4-Isopropyltoluene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
2,2-Dichloropropane	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Benzene	<b>EPA 8260-VOC</b>	D018	Organic Chemistry-VOC	<b>9,16</b>
Bromobenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Bromochloromethane	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Bromodichloromethane	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Bromoform	<b>EPA 8260-VOC</b>	U225	Organic Chemistry-VOC	<b>9</b>
Bromomethane	<b>EPA 8260-VOC</b>	U029	Organic Chemistry-VOC	<b>9</b>
Carbon tetrachloride	<b>EPA 8260-VOC</b>	D019	Organic Chemistry-VOC	<b>9,16</b>
Chlorobenzene	<b>EPA 8260-VOC</b>	D021	Organic Chemistry-VOC	<b>9,16</b>
Chlorodibromomethane	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Chloroethane	<b>EPA 8260-VOC</b>	U045	Organic Chemistry-VOC	<b>9</b>
Chloroform	<b>EPA 8260-VOC</b>	D022	Organic Chemistry-VOC	<b>9,16</b>
Chloromethane	<b>EPA 8260-VOC</b>	U045	Organic Chemistry-VOC	<b>9</b>
cis-1,2-Dichloroethene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
cis-1,3-Dichloropropene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Resource Conservation and Recovery Act (RCRA) – EPA 8260-VOC**

<b>Analyte</b>	<b>Test Code</b>	<b>EPA HW No.</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Dibromomethane	<b>EPA 8260-VOC</b>	U068	Organic Chemistry-VOC	<b>9</b>
Dichlorodifluoromethane	<b>EPA 8260-VOC</b>	U075	Organic Chemistry-VOC	<b>9</b>
Ethylbenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Ethylene dibromide	<b>EPA 8260-VOC</b>	U067	Organic Chemistry-VOC	<b>9</b>
Hexachlorobutadiene	<b>EPA 8260-VOC</b>	D033	Organic Chemistry-VOC	<b>9,16</b>
Isopropylbenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Methyl ethyl ketone*	<b>EPA 8260-VOC</b>	D035	Organic Chemistry-VOC	<b>9,16</b>
Methylene chloride	<b>EPA 8260-VOC</b>	U080	Organic Chemistry-VOC	<b>9</b>
MTBE	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Napthalene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
n-Butylbenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
n-Propylbenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
sec-Butylbenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Styrene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
tert-Butylbenzene	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>
Tetrachloroethene (PCE)	<b>EPA 8260-VOC</b>	D039	Organic Chemistry-VOC	<b>9,16</b>
Toluene	<b>EPA 8260-VOC</b>	U220	Organic Chemistry-VOC	<b>9</b>
trans-1,2-Dichloroethene	<b>EPA 8260-VOC</b>	U079	Organic Chemistry-VOC	<b>9</b>
trans-1,3-Dichloropropene	<b>EPA 8260-VOC</b>	U084	Organic Chemistry-VOC	<b>9</b>
Trichloroethene (TCE)	<b>EPA 8260-VOC</b>	D040	Organic Chemistry-VOC	<b>9,16</b>
Trichlorofluoromethane	<b>EPA 8260-VOC</b>	U121	Organic Chemistry-VOC	<b>9</b>
Vinyl chloride	<b>EPA 8260-VOC</b>	D043	Organic Chemistry-VOC	<b>9,16</b>
Xylene	<b>EPA 8260-VOC</b>	U239	Organic Chemistry-VOC	<b>9</b>
*Add-ons	<b>EPA 8260-VOC</b>		Organic Chemistry-VOC	<b>9</b>

\*Tentative analysis by special request, record specific analyte on test request form and schedule at 801-965-2400



**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Resource Conservation and Recovery Act (RCRA) – EPA 8270-SVOC**

Analyte	Test Code	EPA HW No.	UPHL Unit	Page Number
1,2-Dichlorobenzene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
1,2,4-Trichlorobenzene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
1,3-Dichlorobenzene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
1,4-Dichlorobenzene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
2-Chloronaphthalene	<b>EPA 8270-SVOC</b>	U047	Organic Chemistry-SVOC	<b>10</b>
2-Chlorophenol	<b>EPA 8270-SVOC</b>	U048	Organic Chemistry-SVOC	<b>10</b>
2-Methyl naphthalene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
2-Methyl phenol (o-Cresol)	<b>EPA 8270-SVOC</b>	D023	Organic Chemistry-SVOC	<b>10,16</b>
2-Methyl-4,6-dinitrophenol	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
2-Nitroaniline	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
2-Nitrophenol	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
2,4-Dichlorophenol	<b>EPA 8270-SVOC</b>	U081	Organic Chemistry-SVOC	<b>10</b>
2,4-Dimethylphenol	<b>EPA 8270-SVOC</b>	U101	Organic Chemistry-SVOC	<b>10</b>
2,4-Dinitrophenol	<b>EPA 8270-SVOC</b>	P048	Organic Chemistry-SVOC	<b>10</b>
2,4-Dinitrotoluene	<b>EPA 8270-SVOC</b>	D030	Organic Chemistry-SVOC	<b>10,16</b>
2,4,5-Trichlorophenol	<b>EPA 8270-SVOC</b>	D041	Organic Chemistry-SVOC	<b>10,16</b>
2,4,6-Trichlorophenol	<b>EPA 8270-SVOC</b>	D042	Organic Chemistry-SVOC	<b>10,16</b>
2,6-Dinitrotoluene	<b>EPA 8270-SVOC</b>	U106	Organic Chemistry-SVOC	<b>10</b>
3-Methyl phenol (m-Cresol)	<b>EPA 8270-SVOC</b>	D024	Organic Chemistry-SVOC	<b>10,16</b>
3-Nitroaniline	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
3,3'-Dichlorobenzidine	<b>EPA 8270-SVOC</b>	U073	Organic Chemistry-SVOC	<b>10</b>
4-Bromophenyl phenyl ether	<b>EPA 8270-SVOC</b>	U030	Organic Chemistry-SVOC	<b>10</b>
4-Chloroaniline	<b>EPA 8270-SVOC</b>	P024	Organic Chemistry-SVOC	<b>10</b>
4-Chlorophenyl phenyl ether	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
4-Chloro-3-methyl phenol	<b>EPA 8270-SVOC</b>	U039	Organic Chemistry-SVOC	<b>10</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Resource Conservation and Recovery Act (RCRA) – EPA 8270-SVOC**

Analyte	Test Code	EPA HW No.	UPHL Unit	Page Number
4-Methyl phenol (p-Cresol)	<b>EPA 8270-SVOC</b>	D025	Organic Chemistry-SVOC	<b>10,16</b>
4-Nitroaniline	<b>EPA 8270-SVOC</b>	P077	Organic Chemistry-SVOC	<b>10</b>
4-Nitrophenol	<b>EPA 8270-SVOC</b>	U170	Organic Chemistry-SVOC	<b>10</b>
Acenaphthene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Acenaphthylene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Aniline	<b>EPA 8270-SVOC</b>	U012	Organic Chemistry-SVOC	<b>10</b>
Anthracene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Benzdine	<b>EPA 8270-SVOC</b>	U021	Organic Chemistry-SVOC	<b>10</b>
Benzo (a) anthracene	<b>EPA 8270-SVOC</b>	U018	Organic Chemistry-SVOC	<b>10</b>
Benzo (a) pyrene	<b>EPA 8270-SVOC</b>	U022	Organic Chemistry-SVOC	<b>10</b>
Benzo (b) fluoranthene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Benzo (g,h,i) perylene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Benzo (k) fluoranthene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Benzylbutylphthalate	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Benzyl alcohol	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Benzoic acid	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
bis (2-chloroethyl) ether	<b>EPA 8270-SVOC</b>	U025	Organic Chemistry-SVOC	<b>10</b>
bis (2-chloroethoxy) methane	<b>EPA 8270-SVOC</b>	U046	Organic Chemistry-SVOC	<b>10</b>
bis (2-chloroisopropyl) ether	<b>EPA 8270-SVOC</b>	U027	Organic Chemistry-SVOC	<b>10</b>
bis (2-ethylhexyl) phthalate	<b>EPA 8270-SVOC</b>	U028	Organic Chemistry-SVOC	<b>10</b>
Chrysene	<b>EPA 8270-SVOC</b>	U050	Organic Chemistry-SVOC	<b>10</b>
Dibenzo (a,h) anthracene	<b>EPA 8270-SVOC</b>	U063	Organic Chemistry-SVOC	<b>10</b>
Dibenzofuran	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Diethyl phthalate	<b>EPA 8270-SVOC</b>	U088	Organic Chemistry-SVOC	<b>10</b>
Dimethyl phthalate	<b>EPA 8270-SVOC</b>	U102	Organic Chemistry-SVOC	<b>10</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Resource Conservation and Recovery Act (RCRA) – EPA 8270-SVOC**

Analyte	Test Code	EPA HW No.	UPHL Unit	Page Number
Di-n-butyl phthalate	<b>EPA 8270-SVOC</b>	U069	Organic Chemistry-SVOC	<b>10</b>
Di-n-octyl phthalate	<b>EPA 8270-SVOC</b>	U107	Organic Chemistry-SVOC	<b>10</b>
Fluoranthene	<b>EPA 8270-SVOC</b>	U120	Organic Chemistry-SVOC	<b>10</b>
Fluorene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Hexachlorobenzene	<b>EPA 8270-SVOC</b>	D032	Organic Chemistry-SVOC	<b>10,16</b>
Hexachlorobutadiene	<b>EPA 8270-SVOC</b>	D033	Organic Chemistry-SVOC	<b>10,16</b>
Hexachlorocyclopentadiene	<b>EPA 8270-SVOC</b>	U130	Organic Chemistry-SVOC	<b>10</b>
Hexachloroethane	<b>EPA 8270-SVOC</b>	D034	Organic Chemistry-SVOC	<b>10,16</b>
Ideno (1,2,3-cd) pyrene	<b>EPA 8270-SVOC</b>	U137	Organic Chemistry-SVOC	<b>10</b>
Isophorone	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
n-Nitrosodimethylamine	<b>EPA 8270-SVOC</b>	P082	Organic Chemistry-SVOC	<b>10</b>
n-Nitrosodiphenylamine	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
n-Nitrosodipropylamine	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Naphthalene	<b>EPA 8270-SVOC</b>	U165	Organic Chemistry-SVOC	<b>10</b>
Nitrobenzene	<b>EPA 8270-SVOC</b>	D036	Organic Chemistry-SVOC	<b>10,16</b>
Pentachlorophenol	<b>EPA 8270-SVOC</b>	D037	Organic Chemistry-SVOC	<b>10,16</b>
Phenanthrene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
Phenol	<b>EPA 8270-SVOC</b>	U188	Organic Chemistry-SVOC	<b>10</b>
Pyrene	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>
*Pyridine	<b>EPA 8270-SVOC</b>	D038	Organic Chemistry-SVOC	<b>10,16</b>
*Add-ons	<b>EPA 8270-SVOC</b>		Organic Chemistry-SVOC	<b>10</b>

\*Tentative analysis by special request, record specific analyte on test request form and schedule at 801-965-2400

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Safe Drinking Water Act (SDWA)**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
2,4-D	<b>EPA 515.1</b>	Organic Chemistry-Herbicides	<b>11,12</b>
2,4,5-TP (Silvex)	<b>EPA 515.1</b>	Organic Chemistry-Herbicides	<b>11,12</b>
3-Hydroxycarbofuran	<b>EPA 531.1</b>	Organic Chemistry-Carbamates	<b>11,12</b>
Aldicarb (Temik)	<b>EPA 531.1</b>	Organic Chemistry-Carbamates	<b>11,12</b>
Aldicarb sulfone	<b>EPA 531.1</b>	Organic Chemistry-Carbamates	<b>11,12</b>
Aldicarb sulfoxide	<b>EPA 531.1</b>	Organic Chemistry-Carbamates	<b>11,12</b>
Alkalinity	<b>SM 2320B-ALK</b>	Inorganic Chemistry	<b>30,34</b>
Aluminum	<b>EPA 200.8-AL</b>	Inorganic Chemistry-Metals	<b>20</b>
Ammonia	<b>EPA 350.1-NH3</b>	Inorganic Chemistry-Nutrients	<b>22,24</b>
Antimony	<b>EPA 200.8-SB</b>	Inorganic Chemistry-Metals	<b>20</b>
Arsenic	<b>EPA 200.8-AS</b>	Inorganic Chemistry-Metals	<b>19,20</b>
Barium	<b>EPA 200.8-BA</b>	Inorganic Chemistry-Metals	<b>20</b>
Beryllium	<b>EPA 200.8-BE</b>	Inorganic Chemistry-Metals	<b>20</b>
Boron	<b>EPA 200.7-B</b>	Inorganic Chemistry-Metals	<b>20</b>
Bromate	<b>EPA 300.0-BRO3</b>	Inorganic Chemistry	<b>28</b>
Bromide	<b>EPA 300.1-BRIC</b>	Inorganic Chemistry	<b>28</b>
Bromodichloromethane	<b>EPA 524.2-THM</b>	Organic Chemistry-THM	<b>8</b>
Bromoform	<b>EPA 524.2-THM</b>	Organic Chemistry-THM	<b>8</b>
Cadmium	<b>EPA 200.8-CD</b>	Inorganic Chemistry-Metals	<b>20</b>
Calcium	<b>EPA 200.7-CA</b>	Inorganic Chemistry-Metals	<b>20</b>
Carbaryl (Sevin)	<b>EPA 531.1</b>	Organic Chemistry-Carbamates	<b>11,12</b>
Carbofuran (Furadan)	<b>EPA 531.1</b>	Organic Chemistry-Carbamates	<b>11,12</b>
Chlorate	<b>EPA 300.0-CLO3</b>	Inorganic Chemistry	<b>28</b>
Chloride	<b>EPA 325.2-CL</b>	Inorganic Chemistry	<b>28,35</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Safe Drinking Water Act (SDWA)**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Chloride	<b>EPA 300.0-CLIC</b>	Inorganic Chemistry	<b>28,35</b>
Chlorite	<b>EPA 300.0-CLO2</b>	Inorganic Chemistry	<b>28</b>
Chlorodibromomethane	<b>EPA 524.2-THM</b>	Organic Chemistry-THM	<b>8</b>
Chloroform	<b>EPA 524.2-THM</b>	Organic Chemistry-THM	<b>8</b>
Chromium	<b>EPA 200.8-CR</b>	Inorganic Chemistry-Metals	<b>20</b>
Chromium-IV	<b>EPA 218.7-CR6+</b>	Inorganic Chemistry	<b>35</b>
Cobalt	<b>EPA 200.8-CO</b>	Inorganic Chemistry-Metals	<b>20</b>
Color	<b>EPA 110.2-L-COLR</b>	Inorganic Chemistry	<b>17,27</b>
Conductivity	<b>EPA 120.1-COND</b>	Inorganic Chemistry	<b>26,41</b>
Copper	<b>EPA 200.8-CU</b>	Inorganic Chemistry-Metals	<b>20</b>
Corrosivity	<b>EPA 1110 CORR</b>	Inorganic Chemistry	<b>17,26</b>
<i>Cryptosporidium</i>	<b>EPA 1623</b>	Environmental Microbiology	<b>39</b>
Cyanide	<b>EPA 335.4-CNCL</b>	Inorganic Chemistry	<b>29</b>
Dalapon	<b>EPA 515.1</b>	Organic Chemistry-Herbicides	<b>11,12</b>
Dibromoacetic acid	<b>SM 6251B-HAA</b>	Organic Chemistry-HAA	<b>8</b>
Dichloroacetic acid	<b>SM 6251B-HAA</b>	Organic Chemistry-HAA	<b>8</b>
Dinoseb	<b>EPA 515.1</b>	Organic Chemistry-Herbicides	<b>11,12</b>
<i>E. coli</i>	<b>SM 9223B-Colilert</b>	Environmental Microbiology	<b>38</b>
<i>Enterococcus</i>	<b>Enterolert</b>	Environmental Microbiology	<b>40</b>
Fecal Coliform	<b>SM 9221E-EC</b>	Environmental Microbiology	<b>38</b>
Fluoride	<b>EPA 300.0-FLIC</b>	Inorganic Chemistry	<b>30,35</b>
<i>Giardia</i>	<b>EPA 1623</b>	Environmental Microbiology	<b>33</b>
Gross Alpha	<b>EPA 900.0-ALPG</b>	Radiologic Chemistry	<b>41</b>
Gross Beta	<b>EPA 900.0-BETG</b>	Radiologic Chemistry	<b>41</b>
Gross Gamma	<b>EPA 901.1-GAMG</b>	Radiologic Chemistry	<b>43</b>

# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## PROGRAM SPECIFIC INDEXES

### Safe Drinking Water Act (SDWA)

Analyte	Test Code	UPHL Unit	Page Number
Hardness	<b>HARD</b>	Inorganic Chemistry	<b>27,35</b>
Heterotrophic Plate Count	<b>SM 9215B-HPC</b>	Environmental Microbiology	<b>38</b>
Iron	<b>EPA 200.7-FE</b>	Inorganic Chemistry-Metals	<b>20</b>
Lead	<b>EPA 200.8-PB</b>	Inorganic Chemistry-Metals	<b>20</b>
<i>Legionella</i>	<b>SM 9260J</b>	Environmental Microbiology	<b>40</b>
Magnesium	<b>EPA 200.7-MG</b>	Inorganic Chemistry-Metals	<b>20</b>
Manganese	<b>EPA 200.8-MN</b>	Inorganic Chemistry-Metals	<b>20</b>
Mercury	<b>EPA 245.1-HG</b>	Inorganic Chemistry-Metals	<b>19,20</b>
Methomyl	<b>EPA 531.1</b>	Organic Chemistry-Carbamates	<b>11,12</b>
Molybdenum	<b>EPA 200.8-MO</b>	Inorganic Chemistry-Metals	<b>20</b>
Monobromoacetic acid	<b>SM 6215B-HAA</b>	Organic Chemistry-HAA	<b>8</b>
Monochloroacetic acid	<b>SM 6215B-HAA</b>	Organic Chemistry-HAA	<b>8</b>
Nickel	<b>EPA 200.8-NI</b>	Inorganic Chemistry-Metals	<b>20</b>
Nitrate+Nitrite	<b>EPA 353.2-NO2+NO3</b>	Inorganic Chemistry-Nutrients	<b>23,24</b>
Nitrite (only)	<b>EPA 353.2-NO2</b>	Inorganic Chemistry-Nutrients	<b>23,24</b>
Odor	<b>EPA 140.1-ODOR</b>	Inorganic Chemistry	<b>17,27</b>
Oxamyl (Vydate)	<b>EPA 531.1</b>	Organic Chemistry-Carbamates	<b>11,12</b>
Pentachlorophenol	<b>EPA 515.1</b>	Organic Chemistry-Herbicides	<b>11,12</b>
Perchlorate	<b>EPA 314.0-CLO4</b>	Inorganic Chemistry	<b>29</b>
pH	<b>EPA 150.1-PH</b>	Inorganic Chemistry	<b>25,34</b>
Phosphate	<b>EPA 365.1-TPO4</b>	Inorganic Chemistry-Nutrients	<b>22,24</b>
Picloram	<b>EPA 515.1</b>	Organic Chemistry-Herbicides	<b>11,12</b>
Potassium	<b>EPA 200.7-K</b>	Inorganic Chemistry-Metals	<b>20</b>
Radium-226	<b>EPA 903.1-226R</b>	Radiologic Chemistry	<b>42</b>
Radium-228	<b>EPA 904.0-226R</b>	Radiologic Chemistry	<b>42</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Safe Drinking Water Act (SDWA)  
Safe Drinking Water Act (SDWA) – EPA 524.2-VOC**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Radon	<b>EPA 913.0-RN</b>	Radiologic Chemistry	<b>43</b>
Selenium	<b>EPA 200.8-SE</b>	Inorganic Chemistry-Metals	<b>20</b>
Silica	<b>EPA 370.1-SIO2</b>	Inorganic Chemistry	<b>34</b>
Silver	<b>EPA 200.8-AG</b>	Inorganic Chemistry-Metals	<b>20</b>
Sodium	<b>EPA 200.7-NA</b>	Inorganic Chemistry-Metals	<b>20</b>
Sulfate	<b>EPA 300.0-SO4C</b>	Inorganic Chemistry	<b>23,35</b>
Sulfide	<b>EPA 376.2-SULI</b>	Inorganic Chemistry	<b>30</b>
Thallium	<b>EPA 200.8-TL</b>	Inorganic Chemistry-Metals	<b>35</b>
TOC	<b>SM 5310B-TOC</b>	Inorganic Chemistry	<b>30</b>
Total Coliform	<b>SM 9223B-Colilert</b>	Environmental Microbiology	<b>38</b>
Total Coliform	<b>SM 9222B-mENDO</b>	Environmental Microbiology	<b>38</b>
Total Dissolved Solids (TDS)	<b>EPA 160.1-TDS</b>	Inorganic Chemistry	<b>34,35</b>
Total Suspended Solids (TSS)	<b>EPA 160.2-TSS</b>	Inorganic Chemistry	<b>34,35</b>
Trichloroacetic acid	<b>SM 6251B-HAA</b>	Organic Chemistry-HAA	<b>8</b>
Turbidity	<b>EPA 180.1-TURB</b>	Inorganic Chemistry	<b>32,35</b>
Uranium	<b>EPA 200.8-UUMS</b>	Inorganic Chemistry-Metals	<b>36,41</b>
UV254	<b>SM 5910B-UV254</b>	Inorganic Chemistry	<b>31</b>
Vanadium	<b>EPA 200.8-V</b>	Inorganic Chemistry-Metals	<b>20</b>
Zinc	<b>EPA 200.8-ZN</b>	Inorganic Chemistry-Metals	<b>20</b>
1,1-Dichloroethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,1-Dichloroethene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,1-Dichloropropene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,1,1-Trichloroethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,1,2-Trichloroethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>



# **Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

## **PROGRAM SPECIFIC INDEXES**

### **Safe Drinking Water Act (SDWA) – EPA 524.2-VOC**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
1,1,1,2-Tetrachloroethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,1,2,2-Tetrachloroethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,2-Dibromo-3-chloropropane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,2-Dichlorobenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,2-Dichloroethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,2-Dichloropropane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,2-Dichlorotoluene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,2,3-Trichlorobenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,2,3-Trichloropropane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,2,4-Trichlorobenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,2,4-Trimethylbenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,3-Dichlorobenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,3-Dichloropropane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,3,5-Trimethylbenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,4-Dichlorobenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,4-Dichlorotoluene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
1,4-Isopropyltoluene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
2,2-Dichloropropane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Benzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Bromobenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Bromochloromethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Bromodichloromethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Bromoform	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Bromomethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Carbon tetrachloride	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>

# **Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

## **PROGRAM SPECIFIC INDEXES**

### **Safe Drinking Water Act (SDWA) – EPA 524.2-VOC**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Chlorobenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Chlorodibromomethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Chloroethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Chloroform	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Chloromethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
cis-1,2-Dichloroethene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
cis-1,3-Dichloropropene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Dibromomethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Dichlorodifluoromethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Ethylbenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Ethylene dibromide	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Hexachlorobutadiene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Isopropylbenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Methylene chloride	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
MTBE	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Napthalene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
n-Butylbenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
n-Propylbenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
sec-Butylbenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Styrene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
tert-Butylbenzene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Tetrachloroethene (PCE)	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Toluene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
trans-1,2-Dichloroethene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
trans-1,3-Dichloropropene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES****Safe Drinking Water Act (SDWA) – EPA 524.2-VOC  
Safe Drinking Water Act (SDWA) – EPA 525.2-SVOC**

<b>Analyte</b>	<b>Test Code</b>	<b>UPHL Unit</b>	<b>Page Number</b>
Trichloroethene (TCE)	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Trichlorofluoromethane	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Vinyl chloride	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
Xylene	<b>EPA 524.2-VOC</b>	Organic Chemistry-VOC	<b>9,12</b>
2,4-Dinitrotoluene	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
2,6-Dinitrotoluene	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
4,4'-DDE	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Acetochlor	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Alachlor	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Aldrin	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
alpha-Chlordane	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Atrazine	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
bis (2-ethylhexyl) adipate	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
bis (2-ethylhexyl) phthalate	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Benzo (a) pyrene	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Bromocil	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Butachlor	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Chlorobiphenyl	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Cyanazine	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Dichlorobiphenyl	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Dieldrin	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Endrin	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
EPTC	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
gamma-Chlordane	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>

**Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

**PROGRAM SPECIFIC INDEXES**

**Safe Drinking Water Act (SDWA) – EPA 525.2-SVOC**

Analyte	Test Code	UPHL Unit	Page Number
Heptachlor	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Heptachlor epoxide	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Heptachlorobiphenyl	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Hexachlorobenzene	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Hexachlorobiphenyl	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Hexachlorocyclopentadiene	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Lindane (gamma-BHC)	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Methoxychlor	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Metolachlor	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Metribuzin	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Molinate	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Octachlorobiphenyl	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Pentachlorobiphenyl	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Pentachlorophenol	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Prometon	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Propachlor	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Simazine	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Terbacil	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Tetrachlorobiphenyl	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
trans-Nonachlor	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Trichlorobiphenyl	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Trifluralin	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>
Toxaphene	<b>EPA 525.2-SVOC</b>	Organic Chemistry-SVOC	<b>10,12</b>

**APPENDIX  
INORGANIC AND METAL CHEMISTRY**

Name: **Selenium Hydride**

Test Code: **SM 3114C (modified)**



Grouping: **Metals**

Application: Drinking Water  
Surface Water  
Groundwater

Analyte: Selenium Hydride

Instructions for Collection:

Allow source to flow for a few minutes until water temperature stabilizes. Slowly fill to top of container, be sure not to over fill to prevent loss of preservative.

Required Containers/Volume:

250 mL Metals – plastic bottle

Preservative and Handling:

Add  $\text{HNO}_3$  to pH <2, refrigerate or store on ice and do not allow to freeze, **recommend next day receipt at lab**, holding times, 180 days to analyze

Method Technology:

Hydride AA

## **INSTRUCTIONS**

### **USE OF TEST REQUEST FORMS**

When filling out test request forms the following fields of information must be completed for proper identification of samples at time of receipt:

**System Name / Agency Name**

**System Number / Agency Code**

**Cost Code / Project Code**

**Contact Information**

(phone number is required)

**Person Submitting Samples**

(point of contact for clarifications)

**Billing Information**

(if submitting samples for first time or if updates are needed)

If you are unsure of your system name and number or agency code, please contact the laboratory prior to sample submission at or for drinking water samples contact the Division of Drinking Water at 801-536-4200. If you do not have a system number and have not submitted samples to the laboratory before, a system number will be assigned at time of sample receipt.

Short lists of Agency Codes, Cost Codes / Project Codes, Sample Source Codes, and Sample Type Codes are referenced on page 72. Source and Type codes may be recorded in comments field. If you do not know which cost code or project code applies to your samples, please contact the laboratory.

### **USE OF CHAIN OF CUSTODY**

When submitting chain of custody samples, please complete the following steps to ensure proper preservation of sample integrity:

**Place Seals on Sample Container Lid or Cap or over the cooler.** (must be initialed and dated at time of collection)

**Identify Continuous Sample Possession**

(signatures for dispatch, courier, relinquish, and so forth)  
located at bottom of chain-of-custody form

**Verify Laboratory Receipt**

(obtain copy of form when signed by DLS staff at time of receipt)

# **Chemical and Environmental Laboratory – Client Services Manual**

State of Utah Public Health Laboratory

## **INSTRUCTIONS**

### **AGENCY CODES**

UDAQ	Air Lead High Volume Study	731
UDWQ	Clean Lakes Project	351
UDWQ	Ground Water Oversight	353
UDWQ	Ground Water Permits	352
UDAQ	Utah Division of Air Quality	321
UDDW	Utah Division of Drinking Water	361
UDRC	Utah Division of Radiation Control	342
UDRC	Utah Division of Radiation Control (Mixed Waste)	343
UDSHW	Utah Division of Solid and Hazardous Waste	365/465
UDER	Utah Division of Environmental Response & Remediation	367
UDWQ	Utah Division of Water Quality (General)	350

### **Public Drinking Water Systems, Billing Direct to Customer**

361B

### **Alliance Drinking Water Systems**

366

BRHL	Bear River District Health Department	900
CUHR	Central Utah Public Health Department	900
DCHD6	Davis County Health Department	900
SLWQ	Salt Lake Valley Health Department	900
SEHP	Southeast Utah Public Health Department	900
SWHC	Southwest Utah Public Health Department	900
SCHD1	Summit County Health Department	900
TCHD	Tooele County Health Department	900
UCHD3	Utah County Health Department	900
WCHD3	Wasatch County Health Department	900
WMHD9	Weber-Morgan Health Department	900

### **General Testing, Billing Direct to Customer**

901

### **SAMPLE SOURCE CODES**

1	Spring
2	Well
3	Stream
4	Lake
7	Pipe/Effluent
20	Influent
22	Primary Clarifier
24	Trickling Filter
26	Secondary Clarifier
30	Above UV
31	Below UV
32	Injection Well
33	Bat Point
14	Other

### **SAMPLE TYPE CODES**

4	Grab Sample	27	Below Thermocline
8	8 hr Composite	28	Indicate Depth
9	24 hr Composite	29	Lake Bottom
15	6 hr Composite	30	Sludge
2	Lake Surface	40	Sediment
22	Lake Indicate Depth	50	Soil
23	Lake Above Thermocline	60	Air
24	Lake Indicate Depth	70	Tissue
25	Lake Mid Depth	80	Blood
26	Indicate Depth	90	Urine



## ***Chemical and Environmental Laboratory – Client Services Manual***

State of Utah Public Health Laboratory

### **Sampling Request Instructions**

To help reduce delays in processing your samples, please fill out one test form for each sample site. One sample site may have multiple bottles, but one form cannot have multiple collection sites. Complete the collection information on all of the bottles that are submitted and place completed form in a resalable bag. The Division of Drinking Water has requested some new fields on the forms—Facility ID and Sampling Point ID. If you are not certain what to put in those fields, please contact The Division of Drinking Water. Also, when making a sample delivery, please include a telephone number that is reachable 24-28 hours after sample has been received.

Please make sure you are using our current shipping address to prevent extra courier charges and delays in your shipment as well as to ensure that your environmental samples get delivered to the proper location once they are received at the lab, use one of these two addresses listed below. If the samples need to go to someone's attention at the lab, please contact them prior to shipping so they can coordinate with sample receiving.

Address options for shipping environmental samples:

FedEx or UPS	USPS
Utah Public Health Laboratory Attn BCES Sample Receiving 4431 S 2700 W Taylorsville, UT 84129	Utah Public Health Laboratory Attn BCES Sample Receiving PO Box 144300 Salt Lake City, UT 84131-9988
Utah Public Health Laboratory 4431 S 2700 W Rm 149 Taylorsville, UT 84129	Utah Public Health Laboratory Rm 149 PO Box 144300 Salt Lake City, UT 84131-9988

If there are any questions, call 801-965-2405. Thank you.

David Dick  
Environmental Chemistry Sample Receiving  
Utah Public Health Laboratory  
4431 S 2700 W  
Taylorsville, UT 84129

Our business hours are Monday-Friday from 8 A.M. to 5 P.M.

# Chemical and Environmental Laboratory – Client Services Manual

State of Utah Public Health Laboratory

## INDEX

<u>Analyte</u>	<u>Method Number</u>	<u>Page</u>
Alkalinity	EPA 200.8-AL	30
Ammonia	EPA 350.1-NH3	22
Annual Inorganics and Metals (18 parameters)	Type 9 (Primary Inorganics and Metals Chemistry)	18
Arsenic	EPA 200.8-T-AS, EPA 200.8-D-AS	19
Biological Oxygen Demand (BOD)	EPA 405.1-BOD	32
Bromate	EPA 300.1-BRO3	28
Bromide	EPA 300.1-BRIC	28
BTEX N	SM 8260B BTEX N	13
Carbamates	EPA 531.1	11,12
Chemical Oxygen Demand (COD)	EPA 410.4-COD	33
Chemistry Types (analytical groupings)	Total-Chemistry Type 2, Type 3, Type 9	35
Chlorate	EPA 300.1-CLO3	28
Chloride	EPA 325.2-CL	28, 35
Chlorite	EPA 300.1-CLO2	28
Chlorophyll-A	SM 10200H-CH-A	33
Chromium-VI	EPA 218.7	36
Color	EPA 110.2-L-COLR	17, 27
Conductivity	EPA 120.1-COND	26, 41
Corrosivity	EPA 1110-CORR	17, 26
<i>Cryptosporidium</i> and <i>Giardia</i>	EPA 1623-Protozoa	39
Cyanide	EPA 335.4-CNCL	29,35
<i>Enterococcus</i>	Enterolert	40
Fluoride	4500C-F	30,35
Geosmin and MIB	EPA 525.2-ODOR	44
Glycols	EPG	15
Gross Alpha	EPA 900.0-ALPG	41
Gross Beta	EPA 900.0-BETG	41
Gross Gamma	EPA 901.1-GAMG	43
Haloacetic Acids (HAA)	EPA 6251-HAA	8
Hardness	HARD	27,35
Herbicides	EPA 515.1-Herbicides, EPA 8151-Herbicides	11,12
Heterotrophic Plate Count	SM 9215B-HPC	38
Lead (air samples)	EPA 200.8-T-Pb	44
Lead and Copper	EPA 200.8-Type 8 Metals (T-PB and T-CU)	18
<i>Legionella</i>	SM 9260J	40
Mercury	EPA 245.1-T-HG, EPA 245.1-D-HG, EPA 7471B (T-HG or D-HG)	19,20
Metals, Total and Dissolved (analytical groupings)	Total-Metals Type 7, Type 9, RCRA Types, Filtered-Metals Type 3, Type 4, 6010-Metals (Soil)	20

## **Chemical and Environmental Laboratory – Client Services Manual**

### State of Utah Public Health Laboratory

New Drinking Water Source (46 parameters)	Type 7 (Total Inorganics and Metals Chemistry)	17
Nitrate and Nitrite	EPA 353.2-NO <sub>2</sub> +NO <sub>3</sub>	23,24
Nitrite (only)	353.2-NO <sub>2</sub>	22
Nitrogen, Total Dissolved (DTN)	SM 4500N	25
Nitrogen, Total Kjeldahl (TKN)	EPA 351.4	25
Nutrients, Total and Dissolved (analytical groupings)	Total-Nutrients Type 2, Type 3, Type 4, Type 6; Filtered Nutrients Type 9	24
Odor (Threshold Odor Number, TON)	EPA 140.1-ODOR	17,27
Oil and Grease	EPA 1664-O/G	14
Perchlorate	EPA 314.0-CLO <sub>4</sub>	29
Periphyton	SM 10300C (modified)	33
Pesticides and SVOCs	EPA 525.2-SVOC, EPA 625-SVOC, EPA 8270-SVOC	10
pH	EPA 150.1-PH	25,34
Phase II / Phase V		12
Phosphate	EPA 365.1-TPO <sub>4</sub>	21,23
Polychlorinated Biphenyls and Organochlorine Pesticides	EPA 608-PCB/OcPest, EPA 8081-OcPest, EPA 8082-PCBs	14
Radium-226	EPA 903.1-226R	42
Radium-228	EPA 904.0-228R	42
Radon	EPA 913.0-RN	43
Settable Solids (SS)	EPA 160.5-SS	34
Silica	EPA 370.1-SIO <sub>2</sub>	34
Sulfate	EPA 300.0-SO <sub>4</sub> C or EPA 375.2-SO <sub>4</sub>	23,35
Sulfide	EPA 376.2-SULI	30
Surfactants	SM 5540C-Surfactants	15
Total Coliform, Fecal Coliform and <i>E. coli</i>	SM 9223B-Colilert and SM 9222B-mENDO	38
Total Dissolved Solids (TDS)	EPA 160.1-TDS	34,35
Total Organic Carbon (TOC)	SM 5310B-TOC	31
Total Petroleum Hydrocarbons (TPH)	EPA 8015-TPH	13
Total Suspended Solids (TSS)	EPA 160.2-TSS	34,35
Total Volatile Solids (TVS)	EPA 160.4-TVS	34
Toxic Characteristic Leaching Procedure (TCLP-Metals)	TCLP-Metals	37
Toxic Characteristic Leaching Procedure (TCLP-Organics)	TCLP-Organics	16
Trihalomethanes (THM)	EPA 524.2-THM	8
Turbidity	EPA 180.1-TURB	32,35
UCMR3	EPA 524.3, EPA 522, EPA 200.8, EPA 300.1, EPA 218.7, EPA 537	45,46
Uranium	EPA 200.8-UUMS	36,41
UV254	SM 5910B-UV254	31
VOCs	EPA 524.2-VOC, EPA 624-VOC, EPA 8260-VOC	9

\*indicates tests that require pre-analysis scheduling at 801.965.2508



# Water Bacteriological Analysis Test Request Form

## Utah Public Health Laboratory

4431 S 2700 W Taylorsville, UT 84129-8600

801 965 2400 Fax 801 969 3238

<http://health.utah.gov/lab/chemistry>

Please fill out this form using block letters and with a black or blue pen  
Do not attach this form to the sample.

### LAB USE SECTION

LAB NUMBER

Received Date and Time Stamp

Analyzed Date and Time Stamp

Customer Number	Public Water System	Facility ID	Sampling Point ID
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

System/Agency Name

Collection Point Description / Project Name (if applicable)

Collection Point Description (Continued)

Please fill in the circles or check boxes next to the appropriate option  
DO NOT make any other marks. If you need assistance with this form call 801 965 2405.

- Reporting
- ☐ State Drinking Water Compliance Samples
- ☐ Private Investigative
- ☐ Repeat
- ☐ Downstream Within 5 Of Original Sample
- ☐ Near First Service Connection
- ☐ Original Site
- ☐ Upstream Within 5 Of Original Sample
- ☐ TCR Repeat ONLY
- ☐ TCR Repeat AND GW Source
- ☐ Other

Repeat Continued  
Original Sample Number

Original Sample Collection Date

-   -

M M D D Y Y Y Y

Chlorine Residual ppm

- ☐ Drinking Water
- ☐ Pool/Spa
- ☐ Raw Water
- ☐ Deionized or RO
- ☐ Other

GWR Trigger Source Sample(s)

Facility ID (Source #):

Representative Site: SSG001

Sample #  of  source(s)

- ☐ Sampled all sources in use at time of positive
- ☐ Other source(s) not in use

Name of Wholesaler

Date Wholesaler notified

M M D D Y Y Y Y

Collector's First Name

Collector's Last Name

Collection Date

M M D D Y Y Y Y

Collection Time

AM PM

Collector's Comments

First Name

Last Name

Address

State Zip

Phone Fax

City

State Zip

Phone Fax

Sample Receipt Conditions	<input type="checkbox"/> Health Department Investigative
Temperature	Ice
<input type="text"/>	<input type="radio"/> Yes
	<input type="radio"/> No
Contact Numbers	
USL: Public Health - Environmental Microbiology	(801) 965-2400
State Division of Drinking Water	(801) 536-4200
Contact Your Local Health Department for Pool, Spa, and Hot Tub Information	
Unsatisfactory Sample	
<input type="checkbox"/> Exceeded Holding Time	
<input type="checkbox"/> Out of Date Container	
<input type="checkbox"/> Not State Lab Container	
<input type="checkbox"/> Container is Broken/Leaking	
<input type="checkbox"/> Incomplete Documentation	
<input type="checkbox"/> Other	
Please Submit a New Sample	
Matrix Testing	Additional Testing
<input type="checkbox"/> Well 51	
<input type="checkbox"/> Well 51 and HPC	
<input type="checkbox"/> Well 97	
<input type="checkbox"/> HPC	<input type="checkbox"/> Well 51
LAB USE SECTION	

Pool/ Spa/ Hot tub samples will have a colilert test and a heterotrophic plate count performed on each sample. All other samples will have only a colilert test performed unless specified in "other". The colilert test consists of coliform and *E. coli* analysis.

#### INSTRUCTIONS FOR COLLECTING WATER SAMPLES

1. Do not rinse bottle or touch the lip of bottle.
2. Use only approved containers.
3. Return sample to lab within 24 hours of collection and refrigerate or hold on ice until delivery, also do not allow to freeze. Preferably hold sample at less than 10 degrees Celsius (50 degrees Fahrenheit)
4. Collect sample by removing aerator from tap and letting water run for 2-3 minutes. Fill bottle above the 100 ml line
5. If collecting sample from lake, pond, or type of source water, submerge the bottle, forcing it forward with an even slow motion
6. Select sampling point that will be representative of the system being tested
7. Fill out test request form completely.

#### STATE OF UTAH COLIFORM REGULATIONS (FOR DRINKING WATER ONLY)

For routine sample which are total coliform positive

1. System must collect the number of repeat samples indicated below for each total coliform positive result
 

Population	# of repeat samples
25-1,000	4
>1,000	3
2. The repeat samples must taken within 24 hours of the original positive sample and the repeat test request must indicate the lab number and date of the original positive sample. Specific locations of repeat samples are as follows:
  - a. within 5 service connections upstream
  - b. within 5 service connections downstream
  - c. at the original sample site
3. Additional samples are required for the next month's sampling. The number of additional samples are as follows
 

Population	# of routine	# of additional samples
25-1,000		4
1,000-2,500	2	3
2,500-3,300	3	2
3,300-4,100	4	
>4,100	5 or more	none

For *E. coli* positive samples and repeat samples resulting in total coliform positive

1. If either the original routine sample or any of the repeat samples are fecal coliform positive for *E. coli*, an acute violation has occurred and public notice is required within 72 hours.
2. If both the original routine sample and all repeat samples are total coliform positive, a non-acute violation has occurred and public notice is required within 14 days

Pool/ Spa/ Hot tub samples will have a colilert test and a heterotrophic plate count performed on each sample. All other samples will have only a colilert test performed unless specified in "other". The colilert test consists of coliform and *E. coli* analysis.

#### INSTRUCTIONS FOR COLLECTING WATER SAMPLES

1. Do not rinse bottle or touch the lip of bottle.
2. Use only approved containers.
3. Return sample to lab within 24 hours of collection and refrigerate or hold on ice until delivery, also do not allow to freeze. Preferably hold sample at less than 10 degrees Celsius (50 degrees Fahrenheit)
4. Collect sample by removing aerator from tap and letting water run for 2-3 minutes. Fill bottle above the 100 ml line
5. If collecting sample from lake, pond, or type of source water, submerge the bottle, forcing it forward with an even slow motion
6. Select sampling point that will be representative of the system being tested
7. Fill out test request form completely.

#### STATE OF UTAH COLIFORM REGULATIONS (FOR DRINKING WATER ONLY)

For routine sample which are total coliform positive

1. System must collect the number of repeat samples indicated below for each total coliform positive result
 

Population	# of repeat samples
25-1,000	4
>1,000	3
2. The repeat samples must taken within 24 hours of the original positive sample and the repeat test request must indicate the lab number and date of the original positive sample. Specific locations of repeat samples are as follows:
  - a. within 5 service connections upstream
  - b. within 5 service connections downstream
  - c. at the original sample site
3. Additional samples are required for the next month's sampling. The number of additional samples are as follows
 

Population	# of routine	# of additional samples
25-1,000	1	4
1,000-2,500	2	3
2,500-3,300	3	2
3,300-4,100	4	
>4,100	5 or more	none

For *E. coli* positive samples and repeat samples resulting in total coliform positive

1. If either the original routine sample or any of the repeat samples are fecal coliform positive for *E. coli*, an acute violation has occurred and public notice is required within 72 hours.
2. If both the original routine sample and all repeat samples are total coliform positive, a non-acute violation has occurred and public notice is required within 14 days



# CHAIN OF CUSTODY

## Utah Public Health Laboratory Chemical and Environmental Services

4431 S 2700 W Taylorsville, UT 84129-8600

801 965 2400 Fax 801 969 3238

<http://health.utah.gov/lab/chemistry>

- ☐ Hand Delivered  
☐ Shipped Samples  
☐ Cooler Returned

System/Agency Name:		System/Agency Number:		Cost/Project Code:		REQUESTED TESTS					Received Date and Time:			
REPORTING/CONTACT					BILLING (list if different)					Receipt temperature	Receipt pH	Sample Receipt Conditions		
Attn: _____					Special Code: _____							Yes	No	
Address: _____					Attn: _____							<input type="checkbox"/>	<input type="checkbox"/>	Documentation complete
City, State, Zip: _____					Address: _____							<input type="checkbox"/>	<input type="checkbox"/>	Proper containers and in-date
Phone: _____					City, State, Zip: _____							<input type="checkbox"/>	<input type="checkbox"/>	Containers intact
Fax: _____					Phone: _____							<input type="checkbox"/>	<input type="checkbox"/>	Within holding time
Email: _____					Fax: _____					<input type="checkbox"/>	<input type="checkbox"/>	Coolant		
Submitted By: _____										<input type="checkbox"/>	<input type="checkbox"/>	Temperature within-range		
										<input type="checkbox"/>	<input type="checkbox"/>	Acceptable pH	<input type="checkbox"/> N/A	
										<input type="checkbox"/>	<input type="checkbox"/>	Custody Seals Intact		
COLLECTION POINT DESCRIPTION		Collectors Initials	Collection Date (mm/dd/yy)	Collection Time (24 hr)	COMMENTS								LAB NUMBER	

Dispatched By:	Date and Time:	Courier Company Name:	Invoice/Airbill #:
Relinquished By:	Date and Time:	Received by:	Date and Time:
Relinquished to USL:PH by:	Date and Time:	Received at USL:PH by:	Date and Time: